

Datasheet: HCA182

BATCH NUMBER 147766

Description:	HUMAN ANTI BEVACIZUMAB
Specificity:	BEVACIZUMAB
Other names:	AVASTIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD17976
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in <i>E. coli</i> . The antibody is tagged with a V5 tag and a double extended Strep-tag. This antibody is supplied as a liquid.
Preparation	StrepTactin affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Bevacizumab

Specificity **Human Anti-Bevacizumab Antibody**, clone AbD17976 is a paratope specific, recombinant, anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug bevacizumab and inhibits the binding of the drug to its target, vascular endothelial growth factor A (VEGF-A). The antibody can be used to measure the levels of free bevacizumab and biosimilar products using bioanalytical assays.

A pair of anti-bevacizumab antibodies can be used to develop a pharmacokinetic (PK) bridging ELISA to measure free drug. This antibody, in monovalent Fab format, is recommended as the capture antibody, paired with Anti-Bevacizumab Antibody in full immunoglobulin format, clone AbD16748_hlgG1 ([HCA184P](#)) as the detection antibody.

Bevacizumab (reference product branded as Avastin) is a humanized monoclonal antibody (IgG1/kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis. Bevacizumab is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

[View a summary of all Anti-Bevacizumab Antibodies.](#)

Affinity The monovalent intrinsic affinity of this antibody was measured as $K^D=0.4$ nM by real time, label-free molecular interaction analysis on immobilized bevacizumab.

ELISA This product may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with [HCA184P](#) as the detection reagent.
Protocol: [PK bridging ELISA to measure free drug](#)

References

1. McFarland, T.J. *et al.* (2015) BEVACIZUMAB LEVELS IN BREAST MILK AFTER LONG-TERM INTRAVITREAL INJECTIONS. [Retina. 35 \(8\): 1670-3.](#)
2. Iwamoto, N. *et al.* (2018) Antibody drug quantitation in coexistence with anti-drug antibodies on nSMOL bioanalysis. [Anal Biochem. 540-541: 30-7.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Acknowledgements Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.

Avastin is a trademark of Genentech, Inc.

Health And Safety Information Material Safety Datasheet documentation #10094 available at: <https://www.bio-rad-antibodies.com/SDS/HCA182>
10094

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Biotin](#), [HRP](#)

Recommended Useful Reagents

[HUMAN ANTI BEVACIZUMAB \(HCA184\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA185\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371591:200612'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)