

Datasheet: HCA184

BATCH NUMBER 155646

Description:	HUMAN ANTI BEVACIZUMAB
Specificity:	BEVACIZUMAB
Other names:	AVASTIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD16748_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Bevacizumab
Specificity	Human Anti-Bevacizumab Antibody, clone AbD16748_IgG1 is a paratope-specific,

recombinant anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody drug bevacizumab and inhibits the binding of the drug to its target, vascular endothelial growth factor A (VEGF-A). The antibody can be used to measure the levels of bevacizumab and biosimilars in serum from patients.

A pair of anti-bevacizumab antibodies can be used to develop a pharmacokinetic (PK) bridging ELISA to measure free drug. This antibody is recommended as the detection antibody, paired with Anti-Bevacizumab Antibody in monovalent Fab format, clone AbD17976 ([HCA182](#)) as the capture antibody.

Bevacizumab (reference product branded as Avastin) is a humanized monoclonal antibody (IgG1/kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis. Bevacizumab is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

[View a summary of all Anti-Bevacizumab Antibodies.](#)

Purity	97.8% (capillary electrophoresis, reducing conditions).
---------------	---------------------------------------------------------

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=2.1\text{nM}$ by real time, label-free molecular interaction analysis on immobilized Bevacizumab.
-----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------

ELISA	This product may be used in a direct ELISA or as a detection reagent in a sandwich ELISA, when conjugated to HRP, together with HCA182 as the capture reagent. Protocol: PK bridging ELISA to measure free drug This antibody is fully human and can be used as a reference standard in an ADA assay. Protocol: ADA bridging ELISA
--------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

References	1. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
-------------------	-------------------------------------------------------------------------------------------------------------------------------------

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Avastin is a trademark of Genentech, Inc.
-------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA184
--------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------

10094

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA182\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA185\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M371593:200612'

Printed on 03 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)