

Datasheet: HCA184P

BATCH NUMBER 163160

Description:	HUMAN ANTI BEVACIZUMAB:HRP
Specificity:	BEVACIZUMAB
Other names:	AVASTIN
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD16748_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL? phage display library and expressed in a human cell line.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 1.0 mg/ml
Immunogen	Bevacizumab
Specificity	Human Anti-Bevacizumab Antibody, clone AbD16748_IgG1 is a paratope-specific,

recombinant anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody drug bevacizumab and inhibits the binding of the drug to its target, vascular endothelial growth factor A (VEGF-A). The antibody can be used to measure the levels of bevacizumab and biosimilars in serum from patients.

A pair of anti-bevacizumab antibodies can be used to develop a pharmacokinetic (PK) bridging ELISA to measure free drug. This antibody is recommended as the detection antibody, paired with Anti-Bevacizumab Antibody in monovalent Fab format, clone AbD17976 ([HCA182](#)) as the capture antibody.

Bevacizumab (reference product branded as Avastin) is a humanized monoclonal antibody (IgG1/kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis. Bevacizumab is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

[View a summary of all Anti-Bevacizumab Antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=2.1\text{nM}$ by real time, label-free molecular interaction analysis on immobilized Bevacizumab.
ELISA	This product may be used in a direct ELISA or as a detection reagent in a sandwich ELISA together with HCA182 as the capture reagent. Protocol: PK bridging ELISA to measure free drug
References	1. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Storage	This product is shipped frozen. When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C , if available, otherwise store at -20°C . It is not recommended to keep aliquots at 4°C for more than one week.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA184P 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HUMAN ANTI BEVACIZUMAB \(HCA182\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M404888:220907'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)