

Datasheet: HCA184P

BATCH NUMBER INN1709

Description:	HUMAN ANTI BEVACIZUMAB:HRP
Specificity:	BEVACIZUMAB
Other names:	AVASTIN
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD16748_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 0.1 mg/ml
Immunogen	Bevacizumab
Specificity	Human Anti-Bevacizumab Antibody, clone AbD16748_IgG1 is a paratope-specific,

recombinant anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody drug bevacizumab and inhibits the binding of the drug to its target, vascular endothelial growth factor A (VEGF-A). The antibody can be used to measure the levels of bevacizumab and biosimilars in serum from patients.

A pair of anti-bevacizumab antibodies can be used to develop a pharmacokinetic (PK) bridging ELISA to measure free drug. This antibody is recommended as the detection antibody, paired with Anti-Bevacizumab Antibody in monovalent Fab format, clone AbD17976 ([HCA182](#)) as the capture antibody.

Bevacizumab (reference product branded as Avastin) is a humanized monoclonal antibody (IgG1/kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis. Bevacizumab is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

[View a summary of all Anti-Bevacizumab Antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=2.1\text{nM}$ by real time, label-free molecular interaction analysis on immobilized Bevacizumab.
ELISA	This product may be used in a direct ELISA or as a detection reagent in a sandwich ELISA together with HCA182 as the capture reagent. Protocol: PK bridging ELISA to measure free drug
References	1. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Storage	Store at -70°C . Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA184P 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HUMAN ANTI BEVACIZUMAB \(HCA182\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371595:200612'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)