

Datasheet: HCA189

Description:	RECOMBINANT HUMAN IgA1 KAPPA
Name:	HUMAN IgA1 KAPPA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD18705_hIgA1
Isotype:	IgA1 Kappa
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgA1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgA prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 0.5 mg/ml
Immunogen	Green Fluorescent Protein (GFP).

External Database**Links****UniProt:**[P01876](#) [Related reagents](#)[P01834](#) [Related reagents](#)**Entrez Gene:**[3493](#) IGHA1 [Related reagents](#)[3514](#) IGKC [Related reagents](#)

Specificity

Recombinant Human IgA1 kappa, clone AbD18705_hIgA1 is a recombinant human IgA1 antibody with a kappa light chain.

Recombinant Human IgA1 kappa, clone AbD18705_hIgA1 is specific for green fluorescent protein (GFP) and has no known reactivity with mammalian proteins or other antigens. This product is recommended for use as a standard in assays designed to measure IgA1 levels, or as a control antibody when using other human antibodies of the same isotype and subclass.

Clone AbD18705_hIgA1 forms part of a range of human recombinant immunoglobins containing the kappa light chain. A similar range of reagents containing the lambda light chain is also available from Bio-Rad.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

Sold under license of U.S. Patents 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad.

Health And Safety Information

Material Safety Datasheet documentation #10094 available at: 10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

Licensed Use

For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory

For research purposes only

Technical Advice

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products**Recommended Useful Reagents**

[RECOMBINANT HUMAN IgA1 LAMBDA \(HCA172\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To
find a

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M386665:210528'

Printed on 29 Aug 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)