

Datasheet: HCA193

BATCH NUMBER 155645

Description:	RECOMBINANT HUMAN IgG2 KAPPA
Name:	HUMAN IgG2 KAPPA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD18705_hIgG2
Isotype:	IgG2 Kappa
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG2 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 0.5 mg/ml

Immunogen Green Fluorescent Protein (GFP).

External Database

Links

UniProt:

[P01834](#) [Related reagents](#)

[P01859](#) [Related reagents](#)

Entrez Gene:

[3514](#) IGKC [Related reagents](#)

[3501](#) IGHG2 [Related reagents](#)

Specificity

Recombinant Human IgG2 kappa is a recombinant human IgG2 antibody with a kappa light chain and a heavy chain belonging to the [G2m\(.\)](#), formally known as the G2m(n-) allotype.

Recombinant Human IgG2 kappa is specific for Green Fluorescent Protein and has no known reactivity with mammalian proteins or other antigens. This product is recommended for use as a standard in assays designed to measure human IgG2 kappa levels, or as a control antibody when using other human antibodies of the same isotype and subclass.

Recombinant Human IgG2 kappa forms part of a range of human recombinant immunoglobins containing the kappa light chain. A similar range of reagents containing the lambda light chain is also available from Bio-Rad.

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Acknowledgements

Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.

Health And Safety Information

Material Safety Datasheet documentation #10094 available at:
<https://www.bio-rad-antibodies.com/SDS/HCA193>
10094

Licensed Use

For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory

For research purposes only

Technical Advice

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

RECOMBINANT HUMAN IgG2 LAMBDA (HCA108)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M371606:200612'

Printed on 06 Apr 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)