

Datasheet: HCA194

**BATCH NUMBER 1802**

<b>Description:</b>	RECOMBINANT HUMAN IgG3 KAPPA
<b>Name:</b>	HUMAN IgG3 KAPPA
<b>Format:</b>	Purified
<b>Product Type:</b>	Recombinant Protein
<b>Clone:</b>	AbD18705_hlgG3
<b>Isotype:</b>	IgG3 Kappa
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG3 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	Total protein concentration 0.5 mg/ml

**Immunogen** Green Fluorescent Protein (GFP).

---

**External Database Links**

**UniProt:**

[P01834](#) [Related reagents](#)

[P01860](#) [Related reagents](#)

**Entrez Gene:**

[3514](#) IGKC [Related reagents](#)

[3502](#) IGHG3 [Related reagents](#)

---

**Specificity**

**Recombinant Human IgG3 kappa, clone AbD18705\_hlgG3** is a recombinant human IgG3 antibody with a kappa light chain and a heavy chain belonging to the G3m(b1) allotype, previously designated G3m5. Clone AbD18705\_hlgG3 is specific for Green Fluorescent Protein and has no known reactivity with mammalian proteins or other antigens. This product is recommended for use as a standard in assays designed to measure IgG3 levels, or as a control antibody when using other human antibodies of the same isotype and subclass.

Recombinant Human IgG3 kappa, clone AbD18705\_hlgG3 forms part of a range of human recombinant immunoglobulins containing the kappa light chain. A similar range of reagents containing the lambda light chain is also available from Bio-Rad.

---

**Storage**

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10094 available at:  
10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

---

**Licensed Use**

For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

---

**Regulatory**

For research purposes only

---

**Technical Advice**

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Useful Reagents

[RECOMBINANT HUMAN IgG3 LAMBDA \(HCA178\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M371607:200612'

**Printed on 07 Feb 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)