

Datasheet: HCA197

Description:	HUMAN ANTI HUMAN HbA1c
Specificity:	HbA1c
Other names:	HEMOGLOBIN A1C
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD15783_IgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 0.5 mg/ml

Immunogen The immunogens used in the generation of this antibody were:
A glycosylated peptide, corresponding to N-terminal 9 amino acids of haemoglobin beta chain coupled to the carrier proteins BSA and TRF
Purified HbA1c

External Database Links

UniProt:

[P68871](#) [Related reagents](#)

Entrez Gene:

[3043](#) HBB [Related reagents](#)

Specificity

Human anti Human HbA1c antibody, clone AbD15783_IgG1 recognizes HbA1c, hemoglobin glycosylated at the N-terminal valine of the hemoglobin β chain. It does not bind to non-glycosylated hemoglobin. Besides binding to the normal glycosylated hemoglobin, Human anti Human HbA1c antibody, clone AbD15783_IgG1 also recognizes the glycosylated hemoglobin variants HbC, HbE, HbD and to a lesser degree HbS. Elevated HbF is not detected due to the absence of the hemoglobin β chain.

The measurement of HbA1c levels is used to identify the average plasma glucose concentration. In diabetes mellitus, higher levels of HbA1c have been associated with cardiovascular disease and other diseases. Levels of more than 6.5% HbA1c (48 mmol/mol) are used as a criterion to diagnose diabetes. This antibody is also available in a Fab format ([HCA179](#)).

Affinity

The affinity of this antibody, as a monovalent fab measured against the immunogen, was measured as $K_D=19$ nM by real time, label-free molecular interaction analysis using HbA1c-BSA and HbA1c-TRF conjugates as immobilized antigen.

The affinity of this antibody, as a monovalent fab measured against purified HbA1c, was measured as $K_D=58$ nM by real time, label-free molecular interaction analysis using full length HbA1c as immobilized antigen

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at $2-8^{\circ}\text{C}$ for short term use (up to 4 weeks) and store the remaining aliquots at -20°C .

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Health And Safety Information

Material Safety Datasheet documentation #10094 available at:
10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

Licensed Use For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HUMAN ANTI HUMAN HbA1c \(HCA179\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA180\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA181\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA200\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA217\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA218\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA219\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M416667:230302'

Printed on 02 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)