

Datasheet: HCA200

Description:	HUMAN ANTI HUMAN HbA1c
Specificity:	HbA1c
Other names:	Glycated hemoglobin, Hb1c, Hemoglobin A1c
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD15792
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications	<p>This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.</p> <table> <tr> <th></th><th>Yes</th><th>No</th><th>Not Determined</th><th>Suggested Dilution</th></tr> <tr> <td>ELISA</td><td>■</td><td></td><td></td><td></td></tr> </table> <p>Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.</p>					Yes	No	Not Determined	Suggested Dilution	ELISA	■			
	Yes	No	Not Determined	Suggested Dilution										
ELISA	■													
Target Species	Human													
Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library, expressed in <i>E. coli</i> . The antibody is tagged with a V5 tag and a double extended Strep-tag at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.													
Preparation	StrepTactin affinity chromatography													
Source	<i>E.coli</i>													
Buffer Solution	Phosphate buffered saline													
Preservative Stabilisers	0.01% Thiomersal													
Approx. Protein	Total protein concentration 0.5 mg/ml													

Concentrations

Immunogen	The immunogens used in the generation of this antibody were: A glycosylated peptide, corresponding to N-terminal 9 amino acids of haemoglobin beta chain coupled to carrier proteins Purified HbA1c
------------------	---

External Database Links

UniProt:

[P68871](#)

[Related reagents](#)

Entrez Gene:

[3043](#)

HBB

[Related reagents](#)

Specificity	Human anti Human HbA1c antibody, clone AbD15792 recognizes HbA1c, i.e. hemoglobin when glycosylated at the N-terminal valine of the hemoglobin β chain. It does not bind to non-glycosylated hemoglobin. Besides binding to the normal glycosylated hemoglobin, it also recognizes the glycosylated hemoglobin variants HbD and HbE, but does not recognize the HbS or HbC variants. Elevated HbF cannot be detected due to the absence of the hemoglobin β chain.
--------------------	---

The measurement of HbA1c levels is used to identify the average plasma glucose concentration. In diabetes mellitus, higher levels of HbA1c have been associated with cardiovascular and other diseases. Levels of more than 6.5% HbA1c (48 mmol/mol) are used as a criterion to diagnose diabetes.

Affinity	The affinity of this antibody, as measured against the immunogen is $K_D=7$ nM as measured by real time, label-free molecular interaction analysis using HbA1c-BSA and HbA1c-TRF conjugates as immobilized antigen. The affinity of this antibody, as measured against purified HbA1c is $K_D=54$ nM as measured by real time, label-free molecular interaction analysis using full length HbA1c as immobilized antigen.
-----------------	---

ELISA	This product is suitable for use in direct or indirect ELISA applications
--------------	---

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at $2-8^{\circ}\text{C}$ for short term use (up to 4 weeks) and store the remaining aliquots at -20°C . Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
-------------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA20010094
Licensed Use	For <i>in vitro</i> research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Biotin](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M420901:230706'

Printed on 12 Aug 2023