

Datasheet: HCA203

Description:	HUMAN ANTI ADALIMUMAB
Specificity:	ADALIMUMAB
Other names:	HUMIRA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD18654_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.

Preparation

Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml.

Immunogen

Adalimumab

Specificity

Human Anti-Adalimumab Antibody, clone AbD18654_hlgG1, is a high affinity

anti-idiotypic antibody that specifically recognizes the adalimumab monoclonal antibody and inhibits the binding of adalimumab to its target, TNF α . The antibody can be used to measure the levels of adalimumab and biosimilar products in bioanalytical assays.

Clone AbD18654_hlgG1 can be used to develop and calibrate immune response assays to measure the anti-drug antibody (ADA) response in patient serum.

Adalimumab (branded as Humira) is a fully human monoclonal antibody (IgG1/kappa) approved for the treatment of ankylosing spondylitis, Crohn's disease, psoriasis and several forms of arthritis. The antibody is specific for TNF α , and acts by blocking the activation of TNF receptors by TNF α . This inactivation down-regulates the inflammatory reactions associated with autoimmune diseases.

[View a summary of all Anti-Adalimumab Antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=0.15$ nM by real time, label-free molecular interaction analysis on immobilized Adalimumab.
ELISA	This product may be used in an indirect ELISA or, when conjugated to HRP, as a detection reagent in a bridging ELISA together with HCA202 as the capture reagent. Protocol: PK bridging ELISA to measure free drug This antibody is fully human and can be used as a reference standard in an ADA assay. Protocol: ADA bridging ELISA
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. Humira® is a registered trademark of Abbott Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: 10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M391689:211015'

Printed on 17 Jan 2022