

## Datasheet: HCA207

**BATCH NUMBER 0912R**

<b>Description:</b>	HUMAN ANTI ADALIMUMAB (DRUG/TARGET COMPLEX)
<b>Specificity:</b>	ADALIMUMAB (DRUG/TARGET COMPLEX)
<b>Other names:</b>	HUMIRA®
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD18754_hlgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
<b>Immunogen</b>	Adalimumab

<b>Specificity</b>	<p><b>Human anti Adalimumab, clone AbD18754_hlgG1</b>, is an antibody that specifically recognizes the adalimumab monoclonal antibody while complexed with TNF<math>\alpha</math>. As clone AbD18754_hlgG1 is a complex-specific antibody that recognizes adalimumab bound to its target, TNF<math>\alpha</math>, it does not recognize either free adalimumab or unbound TNF<math>\alpha</math>.</p> <p>Clone AbD18754_hlgG1 can be used in assays to measure adalimumab that has been captured via immobilized TNF<math>\alpha</math>. It can also be used to measure adalimumab/TNF<math>\alpha</math> and biosimilar products using bioanalytical assays.</p> <p>Adalimumab (branded as Humira<sup>®</sup>) is a fully human monoclonal antibody (IgG1/kappa) approved for the treatment of ankylosing spondylitis, Crohn's disease, psoriasis and several forms of arthritis. The antibody is specific for TNF<math>\alpha</math>, and acts by blocking the activation of TNF receptors by TNF<math>\alpha</math>. This inactivation down-regulates the inflammatory reactions associated with autoimmune diseases.</p> <p><a href="#">View a summary of all anti-adalimumab antibodies</a></p>
<b>Affinity</b>	The monovalent intrinsic affinity of this antibody was measured as $K_D=67$ nM by real time, label-free molecular interaction analysis on immobilized adalimumab/TNF $\alpha$ .
<b>ELISA</b>	This product may be used in an indirect ELISA on adalimumab/TNF $\alpha$ complexes or as detection antibody, when conjugated to HRP, for adalimumab bound to TNF $\alpha$ . Protocol: <a href="#">PK antigen capture ELISA to measure bound drug exclusively</a>
<b>References</b>	1. Harth, S. <i>et al.</i> (2019) Generation by phage display and characterization of drug-target complex-specific antibodies for pharmacokinetic analysis of biotherapeutics. <a href="#">MAbs. 11 (1): 178-190.</a>
<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Humira <sup>®</sup> is a registered trademark of Abbott Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10094 available at: <a href="https://www.bio-rad-antibodies.com/SDS/HCA20710094">https://www.bio-rad-antibodies.com/SDS/HCA20710094</a>
<b>Licensed Use</b>	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M371624:200612'

**Printed on 12 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)