

Datasheet: HCA207

BATCH NUMBER 156417

Description:	HUMAN ANTI ADALIMUMAB (DRUG/TARGET COMPLEX)
Specificity:	ADALIMUMAB (DRUG/TARGET COMPLEX)
Other names:	HUMIRA®
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD18754_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Adalimumab

Specificity	<p>Human anti Adalimumab, clone AbD18754_hlgG1, is an antibody that specifically recognizes the adalimumab monoclonal antibody while complexed with TNFα. As clone AbD18754_hlgG1 is a complex-specific antibody that recognizes adalimumab bound to its target, TNFα, it does not recognize either free adalimumab or unbound TNFα.</p> <p>Clone AbD18754_hlgG1 can be used in assays to measure adalimumab that has been captured via immobilized TNFα. It can also be used to measure adalimumab/TNFα and biosimilar products using bioanalytical assays.</p> <p>Adalimumab (branded as Humira[®]) is a fully human monoclonal antibody (IgG1/kappa) approved for the treatment of ankylosing spondylitis, Crohn's disease, psoriasis and several forms of arthritis. The antibody is specific for TNFα, and acts by blocking the activation of TNF receptors by TNFα. This inactivation down-regulates the inflammatory reactions associated with autoimmune diseases.</p> <p>View a summary of all anti-adalimumab antibodies</p>
Affinity	<p>The monovalent intrinsic affinity of this antibody was measured as $K_D=67$ nM by real time, label-free molecular interaction analysis on immobilized adalimumab/TNFα.</p>
ELISA	<p>This product may be used in an indirect ELISA on adalimumab/TNFα complexes or as detection antibody, when conjugated to HRP, for adalimumab bound to TNFα. Protocol: PK antigen capture ELISA to measure bound drug exclusively</p>
References	<p>1. Harth, S. <i>et al.</i> (2019) Generation by phage display and characterization of drug-target complex-specific antibodies for pharmacokinetic analysis of biotherapeutics. MAbs. 11 (1): 178-190.</p>
Storage	<p>Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	<p>12 months from date of despatch</p>
Acknowledgements	<p>Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Humira[®] is a registered trademark of Abbott Laboratories.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA20710094</p>
Licensed Use	<p>For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.</p>

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371624:200612'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)