

Datasheet: HCA210P

BATCH NUMBER 161451

Description:	HUMAN ANTI USTEKINUMAB:HRP
Specificity:	USTEKINUMAB
Other names:	Stelara®
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD17827_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. Conjugated to horseradish peroxidase (HRP) - liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Ustekinumab
Specificity	Human Anti-Ustekinumab Antibody, clone AbD17827_hlgG1 is a paratope specific,

high affinity, anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug ustekinumab and inhibits binding to its target. The antibody can be used to measure the levels of ustekinumab and biosimilar products in bioanalytical assays.

Clone AbD17827_hIgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. Additionally it can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody is recommended as the detection antibody, paired with an anti-ustekinumab antibody in monovalent Fab format as the capture antibody, clone AbD17829 ([HCA208](#)) or clone AbD18334 ([HCA209](#)).

Ustekinumab (Stelara) is a fully human monoclonal antibody drug (IgG1/kappa) that has been approved in US, Canada, and Europe for the treatment of plaque psoriasis and psoriatic arthritis. The FDA approved the medication for the treatment of moderate-severe Crohn's disease in 2016. The drug blocks interleukin IL-12 and IL-23 by binding with high affinity and specificity to the p40 subunit, preventing the binding of these two cytokines to T-lymphocytes, and thereby disrupting the proinflammatory pathway.

[View a summary of all anti-ustekinumab antibodies](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=0.2$ nM by real time, label-free molecular interaction analysis on immobilized ustekinumab.
ELISA	This product may be used in a direct ELISA or as a detection antibody in a bridging ELISA together with AbD17829 (HCA208) or AbD18334 (HCA209) as the capture reagent. Protocol: PK bridging ELISA to measure free drug Protocol: PK bridging ELISA to measure free drug
References	1. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Storage	This product is shipped frozen. When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA210P 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any

form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HUMAN ANTI USTEKINUMAB \(HCA209\)](#)

[HUMAN ANTI USTEKINUMAB \(HCA208\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M395231:220329'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)