

Datasheet: HCA212

**BATCH NUMBER 156411**

<b>Description:</b>	HUMAN ANTI INFLIXIMAB
<b>Specificity:</b>	INFLIXIMAB
<b>Other names:</b>	Remicade
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD17837
<b>Isotype:</b>	HuCAL Fab monovalent
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in <i>E. coli</i> . The antibody is tagged with a V5 tag and a double extended Strep-tag (SAWSHPQFEKGGGSGGGSGGSAWSHPQFEK).
---------------------	--

<b>Preparation</b>	StrepTactin affinity chromatography
--------------------	-------------------------------------

<b>Buffer Solution</b>	Phosphate buffered saline
------------------------	---------------------------

<b>Preservative Stabilisers</b>	0.01% Thiomersal
---------------------------------	------------------

<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
---------------------------------------	----------------------------------

<b>Immunogen</b>	Infliximab
------------------	------------

<b>Specificity</b>	<p><b>Human Anti-Infliximab, clone AbD17837</b>, is an anti-idiotypic antibody that specifically recognizes the infliximab monoclonal antibody and inhibits the binding of infliximab to its target. It can be used in bioanalytical assays to measure the levels of infliximab and biosimilar products, such as Inflectra and Renflexis, in patient samples.</p> <p>Infliximab (branded as Remicade) is a chimeric monoclonal antibody drug (IgG1/kappa) that has been approved for treatment of psoriasis, Crohn's disease, ankylosing spondylitis, psoriatic arthritis, rheumatoid arthritis, and ulcerative colitis.</p> <p>Infliximab is directed against Tumour Necrosis Factor Alpha (TNF<math>\alpha</math>) and acts by blocking the binding of this cytokine to its receptors. Infliximab also induces apoptosis of TNF<math>\alpha</math> expressing T-lymphocytes.</p> <p><a href="#">View a summary of all Anti-Infliximab Antibodies.</a></p>
<b>Affinity</b>	<p>The monovalent intrinsic affinity of this antibody was measured as <math>K_D=3.2</math> nM by real time, label-free molecular interaction analysis on immobilized infliximab.</p>
<b>ELISA</b>	<p>This product may be used in an indirect ELISA or as a capture antibody in an Infliximab bridging ELISA together with HRP-conjugated <a href="#">HCA213</a> or <a href="#">HCA215</a> as the detection reagent.</p> <p>Protocol: <a href="#">PK bridging ELISA to measure free drug</a></p>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	<p>12 months from date of despatch</p>
<b>Acknowledgements</b>	<p>Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.</p> <p>Remicade is a registered trademark of Janssen Biotech Inc. Inflectra is a registered trademark of Hospira, a Pfizer Company. RENFLEXIS is a trademark of Merck Sharp &amp; Dohme Corp.</p>
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10094 available at: <a href="https://www.bio-rad-antibodies.com/SDS/HCA212">https://www.bio-rad-antibodies.com/SDS/HCA212</a> 10094</p>
<b>Licensed Use</b>	<p>For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.</p>
<b>Regulatory</b>	<p>For research purposes only</p>

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Biotin](#), [HRP](#)

### Recommended Useful Reagents

[HUMAN ANTI INFLIXIMAB \(HCA213\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA214\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA215\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA216\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M371632:200612'

**Printed on 12 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)