

Datasheet: HCA212G

BATCH NUMBER 163856

Description:	HUMAN ANTI INFLIXIMAB
Specificity:	INFLIXIMAB
Other names:	Remicade
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD17837
Isotype:	HuCAL Fab monovalent
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library, expressed in <i>E. coli</i> . The antibody is tagged with a V5 tag and a double extended Strep-tag (SAWSHPQFEKGGGSGGGSGGSAWSHPQFEK).
---------------------	---

Preparation	StrepTactin affinity chromatography
--------------------	-------------------------------------

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.01% Thiomersal
---------------------------------	------------------

Approx. Protein Concentrations	Antibody concentration 1.0 mg/ml
---------------------------------------	----------------------------------

Immunogen	Infliximab
------------------	------------

Specificity	<p>Human Anti-Infliximab, clone AbD17837, is an anti-idiotypic antibody that specifically recognizes the infliximab monoclonal antibody and inhibits the binding of infliximab to its target. It can be used in bioanalytical assays to measure the levels of infliximab and biosimilar products, such as Inflectra and Renflexis, in patient samples.</p> <p>Infliximab (branded as Remicade) is a chimeric monoclonal antibody drug (IgG1/kappa) that has been approved for treatment of psoriasis, Crohn's disease, ankylosing spondylitis, psoriatic arthritis, rheumatoid arthritis, and ulcerative colitis.</p> <p>Infliximab is directed against Tumour Necrosis Factor Alpha (TNFα) and acts by blocking the binding of this cytokine to its receptors. Infliximab also induces apoptosis of TNFα expressing T-lymphocytes.</p> <p>View a summary of all Anti-Infliximab Antibodies.</p>
Affinity	<p>The monovalent intrinsic affinity of this antibody is $K_D=3.2$ nM as measured by real time, label-free molecular interaction analysis on immobilized infliximab.</p>
ELISA	<p>This product may be used in an indirect ELISA or as a capture antibody in an Infliximab bridging ELISA together with HRP-conjugated HCA213 or HCA215 as the detection reagent.</p> <p>Protocol: PK bridging ELISA to measure free drug.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	<p>12 months from date of despatch</p>
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>Remicade is a trademark of Janssen Biotech Inc. Inflectra is a trademark of Hospira, a Pfizer Company. RENFLEXIS is a trademark of Merck Sharp & Dohme Corp.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA212G10094</p>
Licensed Use	<p>For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.</p>
Regulatory	<p>For research purposes only</p>

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide STREP-TAG CLASSIC (MCA2489...) [HRP](#)

Mouse Anti Viral V5-TAG (MCA1360...) [Biotin](#), [HRP](#)

Recommended Useful Reagents

[HUMAN ANTI INFLIXIMAB \(HCA213\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA214\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA215\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA216\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M407513:221007'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)