

## Datasheet: HCA213P

<b>Description:</b>	HUMAN ANTI INFLIXIMAB:HRP
<b>Specificity:</b>	INFLIXIMAB
<b>Other names:</b>	REMICADE
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD17841_hlgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/100 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Product Form

Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. Conjugated to horseradish peroxidase (HRP) - liquid.

#### Preparation

Purified IgG prepared by affinity chromatography on Protein A

#### Buffer Solution

Phosphate buffered saline

#### Preservative Stabilisers

0.01% Thiomersal

#### Approx. Protein Concentrations

Total protein concentration 0.1 mg/ml

#### Immunogen

Infliximab

#### Specificity

**Human Anti-Infliximab, clone AbD17841\_hlgG1** is an anti-idiotypic antibody that specifically recognizes the infliximab monoclonal antibody and inhibits the binding of

infliximab to its target. It can be used in bioanalytical assays to measure the levels of infliximab and biosimilar products, such as Inflectra and Renflexis, in patient samples.

Human anti Infliximab, clone AbD17841\_hIgG1 is a fully human antibody and may be used to develop and calibrate immune response assays to measure the anti-drug antibody (ADA) response in patient sera. In the HRP-conjugated format, clone AbD17841\_IgG1 is suitable as a detection antibody in pharmacokinetic (PK) assays.

Infliximab (branded as Remicade) is a chimeric monoclonal antibody drug (IgG1/kappa) that has been approved for treatment of psoriasis, Crohn's disease, ankylosing spondylitis, psoriatic arthritis, rheumatoid arthritis, and ulcerative colitis.

Infliximab is directed against Tumor Necrosis Factor Alpha (TNF $\alpha$ ) and acts by blocking the binding of this cytokine to its receptors. Infliximab also induces apoptosis of TNF $\alpha$  expressing T-lymphocytes.

[View a summary of all Anti-Infliximab Antibodies.](#)

---

<b>Affinity</b>	The monovalent intrinsic affinity of this antibody in an unconjugated format was measured as $K_D=1.8$ nM by real time, label-free molecular interaction analysis on immobilized infliximab.
<b>ELISA</b>	This product may be used in as a detection antibody in an Infliximab bridging ELISA together with <a href="#">HCA212</a> as the capture reagent. Protocol: <a href="#">PK bridging ELISA to measure free drug</a>
<b>References</b>	1. Hernández-Brejido, B. <i>et al.</i> (2016) Standardization of the homogeneous mobility shift assay protocol for evaluation of anti-infliximab antibodies. Application of the method to Crohn's disease patients treated with infliximab. <a href="#">Biochem Pharmacol. Sep 21. pii: S0006-2952(16)30302-1. [Epub ahead of print]</a>
<b>Storage</b>	This product is shipped frozen. When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <a href="http://bio-rad.com/en-us/trademarks">bio-rad.com/en-us/trademarks</a> for details. Remicade® is a registered trademark of Janssen Biotech Inc. Inflectra® is a registered trademark of Hospira, a Pfizer Company. RENFLEXIS™ is a trademark of Merck Sharp & Dohme Corp.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10094 available at: 10094: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf</a>

---

**Licensed Use** For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Useful Reagents

[HUMAN ANTI INFLIXIMAB \(HCA212\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA214\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA215\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA216\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M391708:211015'

**Printed on 19 May 2022**

---

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)