

## Datasheet: HCA216P

**BATCH NUMBER 166161**

<b>Description:</b>	HUMAN ANTI INFLIXIMAB:HRP
<b>Specificity:</b>	INFLIXIMAB
<b>Other names:</b>	REMICADE®
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD19376_hlgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody selected from the HuCAL phage display library and expressed in a human cell line. Conjugated to horseradish peroxidase (HRP) - liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A.
<b>Source</b>	HKB-11
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.01% Thiomersal.
<b>Approx. Protein Concentrations</b>	Total protein concentration 1.0 mg/ml.
<b>Immunogen</b>	Infliximab

**Specificity** **Human anti Infliximab, clone AbD19376\_hlgG1**, is an anti-idiotypic antibody that binds to both free infliximab and to infliximab bound to its target, TNF $\alpha$ . Clone AbD19376\_hlgG1 does not inhibit the binding of infliximab to TNF $\alpha$ . It can be used in bioanalytical assays to measure the levels of infliximab and biosimilar products, such as Inflectra and Renflexis, in patient samples.

Clone AbD19376\_hlgG1 binds to an epitope outside the infliximab binding site. This epitope is shared with other chimeric antibodies that are highly homologous to infliximab including the chimeric antibody cetuximab but not rituximab.

Human anti Infliximab, clone AbD19376\_hlgG1 is a fully human antibody and may be used to develop and calibrate immune response assays to measure the anti-drug antibody (ADA) response in patient sera. In the HRP-conjugated format, clone AbD19376\_IgG1 is suitable as a detection antibody in pharmacokinetic (PK) assays to measure total infliximab levels (free drug, partially bound and fully bound).

Infliximab (branded as Remicade®) is a chimeric monoclonal antibody drug (IgG1/kappa) that has been approved for treatment of psoriasis, Crohn's disease, ankylosing spondylitis, psoriatic arthritis, rheumatoid arthritis, and ulcerative colitis.

Infliximab is directed against Tumor Necrosis Factor Alpha (TNF $\alpha$ ) and acts by blocking the binding of this cytokine to its receptors. Infliximab also induces apoptosis of TNF $\alpha$  expressing T-lymphocytes.

[View a summary of all anti-infliximab antibodies.](#)

---

**Affinity** The monovalent intrinsic affinity of this antibody in an unconjugated format was measured as  $K_D=0.7$  nM by real time, label-free molecular interaction analysis on immobilized infliximab.

---

**ELISA** This product may be used as a detection reagent in a bridging ELISA for PK assay development together with [HCA214](#) or [HCA212](#) as the capture reagent. Protocol: [PK bridging ELISA to measure total drug.](#)

---

**References** 1. Kashiwagi, N. *et al.* (0) Method for measuring anti-drug antibody [US Patent Application US2017315118A1](#)

---

**Storage** This product is shipped frozen.  
When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.

---

**Guarantee** 12 months from date of despatch.

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.  
Remicade is a registered trademark of Janssen Biotech Inc. Inflectra is a registered

trademark of Hospira, a Pfizer Company. RENFLEXIS is a trademark of Merck Sharp & Dohme Corp.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10094 available at:  
<https://www.bio-rad-antibodies.com/SDS/HCA216P>  
10094

---

**Licensed Use**

For *in vitro*. research purposes and for commercial applications for the provision of *in vitro*. testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Technical Advice**

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Useful Reagents

[HUMAN ANTI INFLIXIMAB \(HCA212\)](#)

[HUMAN ANTI INFLIXIMAB \(HCA214\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[RECOMBINANT HUMAN TNF ALPHA \(PHP051\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M421895:230807'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)