

## Datasheet: HCA217

<b>Description:</b>	HUMAN ANTI HUMAN HbA1c
<b>Specificity:</b>	HbA1c
<b>Other names:</b>	HEMOGLOBIN A1C
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD15785_IgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	Total protein concentration 0.5 mg/ml
<b>Immunogen</b>	The immunogens used in the generation of this antibody were: A glycosylated peptide, corresponding to N-terminal 9 amino acids of haemoglobin beta chain

coupled to the carrier proteins BSA and TRF  
Purified HbA1c

---

**External Database  
Links**

**UniProt:**

[P68871](#)    [Related reagents](#)

**Entrez Gene:**

[3043](#)    HBB    [Related reagents](#)

---

**Specificity**

**Human anti Human HbA1c antibody, clone AbD15785\_IgG1** recognizes HbA1c, hemoglobin glycosylated at the N-terminal valine of the hemoglobin  $\beta$  chain. It does not bind to non-glycosylated hemoglobin. Besides binding to the normal glycosylated hemoglobin. Human anti Human HbA1c antibody, clone AbD15785\_IgG1 also recognizes the glycosylated hemoglobin variants HbC, HbE, HbD and to a lesser degree HbS. Elevated HbF is not detected due to the absence of the hemoglobin  $\beta$  chain.

The measurement of HbA1c levels is used to identify the average plasma glucose concentration. In diabetes mellitus, higher levels of HbA1c have been associated with cardiovascular disease and other diseases. Levels of more than 6.5% HbA1c (48 mmol/mol) are used as a criterion to diagnose diabetes. This antibody is also available in a Fab format ([HCA180](#)).

---

**Affinity**

The affinity of this antibody, as a monovalent fab measured against the immunogen, was measured as  $K_D=28$  nM by real time, label-free molecular interaction analysis using HbA1c-BSA and HbA1c-TRF conjugates as immobilized antigen.

The affinity of this antibody, as a monovalent fab measured against purified HbA1c, was measured as  $K_D=65$  nM by real time, label-free molecular interaction analysis using full length HbA1c as immobilized antigen

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at  $-20^{\circ}\text{C}$  on receipt. When thawed, aliquot the sample as needed. Keep aliquots at  $2-8^{\circ}\text{C}$  for short term use (up to 4 weeks) and store the remaining aliquots at  $-20^{\circ}\text{C}$ .

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10094 available at:  
10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

---

**Licensed Use**

For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Useful Reagents

[HUMAN ANTI HUMAN HbA1c \(HCA179\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA180\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA181\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA197\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA200\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA218\)](#)

[HUMAN ANTI HUMAN HbA1c \(HCA219\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M391714:211015'

**Printed on 15 Oct 2021**

---

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)