

Datasheet: HCA228

BATCH NUMBER 160605

Description:	HUMAN ANTI CETUXIMAB
Specificity:	CETUXIMAB
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD19376_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Chimeric monoclonal antibody
Specificity	Human Anti-Cetuximab Antibody, clone AbD19376_hlgG1 , is a non-inhibitory anti-idiotypic antibody that binds to cetuximab. This antibody recognizes an epitope outside the cetuximab binding site and so does not inhibit the binding of the drug to its

target, epidermal growth factor receptor (EGFR). This epitope is shared with other chimeric antibodies that are highly homologous to cetuximab including infliximab but excluding rituximab. The antibody is suitable for the development of bioanalytical assays for cetuximab and biosimilars.

Clone AbD19376_hlgG1 is a fully human recombinant monoclonal with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. In addition, it is suitable as a detection antibody in pharmacokinetic (PK) bridging ELISA, with Anti-Cetuximab Antibody [HCA220](#) as the capture antibody. This combination enables measurement of free and partially bound cetuximab, which is equivalent to total drug when trough serum levels are above 4 µg/ml.

Cetuximab (brand name Erbitux) is a chimeric monoclonal antibody drug (IgG1/kappa) approved for treatment of KRAS wild-type metastatic colorectal cancers and squamous cell carcinoma of the head and neck. This therapeutic antibody is directed against EGFR and inhibits tumor cell proliferation by blocking the interaction of epidermal growth factor (EGF) with its receptor.

[View a summary of all anti cetuximab antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=39$ nM by real time, label-free molecular interaction analysis on immobilized cetuximab.
ELISA	This product may be used in a direct ELISA or when conjugated to HRP, as a detection antibody in a cetuximab bridging ELISA for PK assay development together with HCA220 as the capture reagent. Protocol: PK bridging ELISA to measure total drug
References	<ol style="list-style-type: none">1. Bose, M.E. <i>et al.</i> (2016) Simulated Respiratory Secretion for Use in the Development of Influenza Diagnostic Assays. PLoS One. 11 (11): e0166800.2. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Further Reading	<ol style="list-style-type: none">1. Gomez, D. <i>et al.</i> (2013) Cetuximab therapy in the treatment of metastatic colorectal cancer: the future frontier? Int J Surg. 11 (7): 507-13.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Health And Safety Information Material Safety Datasheet documentation #10094 available at:
<https://www.bio-rad-antibodies.com/SDS/HCA228>
10094

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HUMAN ANTI CETUXIMAB \(HCA220\)](#)

[HUMAN ANTI CETUXIMAB \(HCA221\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M391722:211015'

Printed on 06 May 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)