

Datasheet: HCA228

Description:	HUMAN ANTI CETUXIMAB
Specificity:	CETUXIMAB
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD19376_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Chimeric monoclonal antibody
Specificity	Human Anti-Cetuximab Antibody, clone AbD19376_hlgG1 , is a non-inhibitory anti-idiotypic antibody that binds to cetuximab. This antibody recognizes an epitope outside the cetuximab binding site and so does not inhibit the binding of the drug to its target, epidermal growth factor receptor (EGFR). This epitope is shared with other

chimeric antibodies that are highly homologous to cetuximab including infliximab but excluding rituximab. The antibody is suitable for the development of bioanalytical assays for cetuximab and biosimilars.

Clone AbD19376_hIgG1 is a fully human recombinant monoclonal with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. In addition, it is suitable as a detection antibody in pharmacokinetic (PK) bridging ELISA, with Anti-Cetuximab Antibody [HCA220](#) as the capture antibody. This combination enables measurement of free and partially bound cetuximab, which is equivalent to total drug when trough serum levels are above 4 µg/ml.

Cetuximab (brand name Erbitux) is a chimeric monoclonal antibody drug (IgG1/kappa) approved for treatment of KRAS wild-type metastatic colorectal cancers and squamous cell carcinoma of the head and neck. This therapeutic antibody is directed against EGFR and inhibits tumor cell proliferation by blocking the interaction of epidermal growth factor (EGF) with its receptor.

[View a summary of all anti cetuximab antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=39$ nM by real time, label-free molecular interaction analysis on immobilized cetuximab.
ELISA	This product may be used in a direct ELISA or when conjugated to HRP, as a detection antibody in a cetuximab bridging ELISA for PK assay development together with HCA220 as the capture reagent. Protocol: PK bridging ELISA to measure total drug
References	1. Bose, M.E. <i>et al.</i> (2016) Simulated Respiratory Secretion for Use in the Development of Influenza Diagnostic Assays. PLoS One. 11 (11): e0166800. 2. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Further Reading	1. Gomez, D. <i>et al.</i> (2013) Cetuximab therapy in the treatment of metastatic colorectal cancer: the future frontier? Int J Surg. 11 (7): 507-13.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
Health And Safety	Material Safety Datasheet documentation #10094 available at:

Information	10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HUMAN ANTI CETUXIMAB \(HCA220\)](#)

[HUMAN ANTI CETUXIMAB \(HCA221\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M371650:200612'

Printed on 12 Feb 2021