

Datasheet: HCA237

**BATCH NUMBER 159489**

<b>Description:</b>	HUMAN ANTI OMALIZUMAB (DRUG/TARGET COMPLEX)
<b>Specificity:</b>	OMALIZUMAB DRUG/TARGET COMPLEX
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD20760_hlgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Omalizumab

**Specificity** **Human anti omalizumab, clone AbD20760\_hlgG1** binds specifically to the omalizumab/IgE (drug/target) complex. The antibody does not recognize free omalizumab or free human IgE. It can be used in bioanalytical assays to detect bound omalizumab or

biosimilar products exclusively.

Clone AbD20760\_hlgG1 can be used to develop a pharmacokinetic (PK) antigen capture assay to measure omalizumab captured via immobilized human IgE or to measure the level of drug/IgE complex in preclinical or clinical samples. This product is in a full immunoglobulin format and is available directly labeled with HRP. The antibody is also available in unlabeled Fab format ([HCA238](#)).

Omalizumab (brand name Xolair<sup>®</sup>) is a recombinant DNA-derived humanized IgG1 kappa monoclonal antibody used in the treatment of patients with moderate or severe asthma who have demonstrated a positive allergy skin test and whose symptoms are not controlled by inhaled corticosteroids.

Allergic asthma is mediated by IgE released by B cells in response to allergen. Circulating IgE binds to the high-affinity IgE Fc receptor (FcεRI) expressed on basophils and mast cells, triggering the release of histamine, leukotrienes and other mediators associated with the pathophysiology of asthma. Omalizumab is directed against the Fc region of human immunoglobulin E (IgE). By binding to circulating IgE at the site of FcεRI binding, this therapeutic antibody prevents the interaction of IgE with its receptor thus limiting mediator release. Treatment with omalizumab has also been demonstrated to reduce the expression of FcεRI on mast cells and basophils, providing additional clinical benefit.

[View a summary of all anti-omalizumab antibodies](#)

---

<b>Affinity</b>	The monovalent intrinsic affinity of this antibody was measured as $K_D = 0.58$ nM by real time, label free molecular interaction analysis on immobilized omalizumab in complex with human IgE.
<b>ELISA</b>	This product may be used in a direct ELISA to detect the omalizumab/IgE complex or to detect omalizumab bound to immobilized IgE. Protocol: <a href="#">PK antigen capture ELISA to measure bound drug exclusively</a>
<b>References</b>	1. Harth, S. <i>et al.</i> (2019) Generation by phage display and characterization of drug-target complex-specific antibodies for pharmacokinetic analysis of biotherapeutics. <a href="#">MAbs. 11 (1): 178-190.</a>
<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.  Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	Sold under license of U.S. Patents 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad. Xolair <sup>®</sup> is a registered trademark of Novartis.

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10094 available at: <a href="https://www.bio-rad-antibodies.com/SDS/HCA237">https://www.bio-rad-antibodies.com/SDS/HCA237</a> 10094
<b>Licensed Use</b>	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a>

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)  
[RECOMBINANT HUMAN IgE KAPPA \(HCA190\)](#)  
[RECOMBINANT HUMAN IgE LAMBDA \(HCA171\)](#)  
[HUMAN ANTI OMALIZUMAB \(DRUG/TARGET COMPLEX\) \(HCA238\)](#)  
[HUMAN ANTI OMALIZUMAB \(HCA236\)](#)  
[HUMAN ANTI OMALIZUMAB \(HCA235\)](#)  
[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M386727:210528'

Printed on 12 Sep 2024