

Datasheet: HCA245

BATCH NUMBER 170671

Description:	HUMAN ANTI GOLIMUMAB (DRUG/TARGET COMPLEX)
Specificity:	GOLIMUMAB DRUG/TARGET COMPLEX
Other names:	SIMPONI
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD20893_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid. Simponi is a trademark of Janssen Biotech, Inc. in the USA and Merck & Co, Inc. in Europe.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml

Immunogen	Golimumab/TNF α complex
Specificity	<p>Human anti golimumab (drug/target complex) antibody, clone AbD20893_hlgG1 is a Type 3 anti-idiotypic antibody that specifically recognizes golimumab when bound to tumor necrosis factor alpha (TNFα) i.e. the drug/target complex. Clone AbD20893_hlgG1 does not recognize unbound golimumab nor free TNFα. Human anti golimumab (drug/target complex) antibody, clone AbD20893_hlgG1 may be used to measure golimumab when captured by immobilized human TNFα or golimumab/TNFα levels in serum from patients.</p> <p>Golimumab, marketed under the brand name Simponi, is a recombinant human IgG1/kappa monoclonal antibody approved for the treatment of rheumatoid arthritis, psoriatic arthritis and ankylosing spondylitis. This therapeutic antibody, directed against TNFα, acts by blocking the binding of TNFα to its receptors, resulting in a down-regulation of the inflammatory response associated with autoimmune diseases.</p> <p>View a summary of all anti-golimumab antibodies</p>
Affinity	The intrinsic affinity of the monovalent form of Human anti Golimumab antibody, clone AbD20893_hlgG1 is $K_D = 53$ nM as measured by real time, label free molecular interaction analysis on immobilized Golimumab in complex with human TNF α .
ELISA	When conjugated to HRP, this product may be used in a direct ELISA to detect the golimumab/TNF α complex or to detect golimumab bound to immobilized TNF α . Protocol: PK antigen capture ELISA to measure bound drug exclusively .
References	1. Harth, S. <i>et al.</i> (2019) Generation by phage display and characterization of drug-target complex-specific antibodies for pharmacokinetic analysis of biotherapeutics. MAbs. 11 (1): 178-90.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA245
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any

form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK003P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK004P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK005P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK006P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M430226:240503'

Printed on 10 Jul 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)