

Datasheet: HCA252G

BATCH NUMBER 161062

Description:	HUMAN ANTI TOCILIZUMAB
Specificity:	TOCILIZUMAB
Other names:	ACTEMRA®, ROACTEMRA®
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD21362
Isotype:	HuCAL Fab monovalent
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA (1)	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) NB. The Rat anti DYKDDDDK-tag antibody, clone 6F7 is not suitable for use with this antibody as clone 6F7 exhibits cross reactivity with the human IL-6R

Product Form

A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL® phage display library, expressed in *E. coli*. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 1.0 mg/ml

Immunogen

Tocilizumab

Specificity	<p>Human Anti-Tocilizumab Antibody, clone AbD21362 is a paratope specific, anti-idiotypic antibody that binds to free tocilizumab but not to human Interleukin-6 Receptor (IL-6R) or the drug/IL-6R complex. The antibody can be used to measure the levels of tocilizumab and biosimilar products in bioanalytical and patient monitoring assays.</p> <p>A pair of anti-tocilizumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody, in monovalent Fab format, is recommended as a capture antibody, paired with an HRP conjugated Anti-Tocilizumab Antibody in full immunoglobulin format, clone AbD21338_hIgG1 (HCA253) as the detection antibody.</p> <p>Tocilizumab (RoActemra/Actemra) is a humanized IgG1/kappa monoclonal antibody approved for the treatment of moderate to severe rheumatoid arthritis and for the treatment of systemic juvenile idiopathic arthritis. The therapeutic antibody is directed against IL-6R and acts by blocking the binding of IL-6R (both the soluble form, sIL-6R, and the membrane bound form, mIL-6R) to IL-6, thus impeding the pro-inflammatory properties of this cytokine.</p> <p>View a summary of all anti-tocilizumab antibodies</p>
Affinity	<p>The monovalent intrinsic affinity of this antibody was measured as $K_D=0.9$ nM by real time, label-free molecular interaction analysis on immobilized tocilizumab.</p>
ELISA	<p>This product may be used in a direct or indirect ELISA or as a capture reagent in a PK bridging ELISA together with HCA253 as the detection reagent.</p> <p>Protocol: PK bridging ELISA to measure free drug</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	<p>12 months from date of despatch</p>
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>Actemra® and RoActemra® are registered trademarks of Hoffman-La Roche and Chugai Pharmaceuticals co., Ltd.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA252G</p> <p>10094</p>

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI TOCILIZUMAB \(HCA253\)](#)

[HUMAN ANTI TOCILIZUMAB \(HCA254\)](#)

[HUMAN ANTI TOCILIZUMAB \(HCA255\)](#)

[HUMAN ANTI TOCILIZUMAB \(HCA256\)](#)

[HUMAN ANTI TOCILIZUMAB \(HCA257\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391753:211015'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)