

Datasheet: HCA262

Description:	HUMAN ANTI PALIVIZUMAB
Specificity:	PALIVIZUMAB
Other names:	SYNAGIS®
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD23967_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody selected from the HuCAL® phage display library and expressed in a human cell line.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Palivizumab

Specificity

Human anti Palivizumab, clone AbD23967_hIgG1, is a paratope specific, high affinity, anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug palivizumab and inhibits it binding to its target. The antibody can be used to measure the levels of palivizumab and biosimilar products in bioanalytical assays.

Clone AbD23967_hIgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a reference standard in an anti-drug antibody (ADA) assay. Additionally the antibody can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody

is recommended as the detection antibody and is available unlabeled or directly conjugated to HRP. Use in a pair with capture antibody clone AbD23967 ([HCA261](#)) in monovalent Fab format.

Palivizumab (Synagis®) is a humanized mouse monoclonal antibody (IgG1/kappa) specific to an epitope of the F protein of respiratory syncytial virus (RSV). It binds the antigenic site A, a highly conserved region on the RSV F protein between amino acids 258 and 275. RSV is known to cause infections of the lower respiratory tract in a clear seasonal pattern during the winter months. Administration of the drug is recommended for premature infants (under 35 weeks gestation), infants with certain congenital heart defects (CHD) or bronchopulmonary dysplasia (BPD) throughout the duration of the RSV season. It is approved by the FDA and EMA as a treatment for respiratory tract infection caused by RSV.

[View a summary of all anti-palivizumab antibodies](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D = 1.5$ nM by real time, label free molecular interaction analysis on immobilized palivizumab.
ELISA	This product may be used in direct or indirect ELISA and in an ELISA bridging assay as a capture or detection antibody together with the same antibody or HCA261 as the paired reagent. It is also suitable as a fully human control or calibrator antibody in an ADA bridging ELISA. Protocol: PK bridging ELISA to measure free drug
References	1. Kashiwagi, N. <i>et al.</i> (2017) Method for measuring anti-drug antibody US Patent Application US20170315118A1
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: 10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HUMAN ANTI PALIVIZUMAB \(HCA261\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M371694:200612'

Printed on 11 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)