

Datasheet: HCA267

BATCH NUMBER 167147

Description:	HUMAN ANTI BOVINE INTERLEUKIN-2
Specificity:	IL-2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD21766
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	
ELISpot	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. The use of Leucoperm (Product Code [BUF09](#)) is recommended for this purpose.

Target Species	Bovine
Product Form	A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL phage display library, expressed in <i>E. coli</i> . The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml.
Immunogen	Fc-fusion protein containing the amino acid sequence 21-155 from bovine interleukin-2.
External Database Links	<p>UniProt: P05016 Related reagents</p> <p>Entrez Gene: 280822 IL2 Related reagents</p>
Synonyms	IL-2
Specificity	<p>Human anti Bovine Interleukin-2, clone AbD21766 recognizes bovine interleukin-2 (IL-2), also known as T-cell growth factor (TCGF). Bovine IL-2 is a 155 amino acid secreted cytokine produced by T-cells following antigenic or mitogenic stimulation. It plays a crucial role in T-cell proliferation and immune regulation.</p> <p>Human anti bovine interleukin-2, clone AbD21766, has been successfully used for the detection and measurement of native and recombinant bovine IL-2 using intracellular flow cytometry, ELISPOT and ELISA.</p> <p>In ELISA and ELISPOT clone AbD21766 functions as the capture antibody in a sensitive assay for bovine IL-2 determination.</p>
Affinity	The intrinsic affinity of the monovalent form of this antibody is K _D =0.3 nM as measured by real time, label-free molecular interaction analysis on recombinant bovine interleukin-2.
Flow Cytometry	Use 10µl of the suggested working dilution to label 1x10 ⁶ cells in 100µl.
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with HCA268P as the detection reagent.
Further Reading	<ol style="list-style-type: none"> 1. Rhodes, S.G. <i>et al.</i> (2014) Use of antigen-specific interleukin-2 to differentiate between cattle vaccinated with Mycobacterium bovis BCG and cattle infected with M. bovis. Clin Vaccine Immunol. 21:39-45. 2. Whelan, A.O., <i>et al.</i> (2011) Development of an Antibody to Bovine IL-2 Reveals Multifunctional CD4 T(EM) Cells in Cattle Naturally Infected with Bovine Tuberculosis PLoS One. 6(12):e29194
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in</p>

frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. His-tag is a trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA267 10040
Licensed Use	For <i>in vitro</i> research purposes only, unless otherwise specified in writing by Bio-Rad
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Negative Controls

[HuCAL Fab-FH NEGATIVE CONTROL \(HCA045\)](#)

Recommended Useful Reagents

[HUMAN ANTI BOVINE INTERLEUKIN-2 \(HCA268\)](#)

[5x ELISA COATING BUFFER \(BUF030A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M421608:230706'

Printed on 01 Mar 2024