

## Datasheet: HCA276

**BATCH NUMBER 166738**

<b>Description:</b>	HUMAN ANTI ETANERCEPT
<b>Specificity:</b>	ETANERCEPT
<b>Other names:</b>	Enbrel®
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD25939
<b>Isotype:</b>	HuCAL Fab monovalent
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library, expressed in *E. coli*. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

### Preparation

Metal chelate affinity chromatography.

### Source

*E.coli*

### Buffer Solution

Phosphate buffered saline.

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

### Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml.

**Specificity**

**Human anti etanercept, clone AbD25939** is an antibody that specifically recognizes etanercept as well as the extracellular ligand-binding portion of the 75 kDa (p75) human tumor necrosis factor receptor 2 (TNFR2). Clone AbD25939 does not inhibit the binding of etanercept to its target tumor necrosis factor alpha (TNF $\alpha$ ). The antibody binds to an epitope outside the TNF $\alpha$  binding site and can be used to measure the levels of etanercept and biosimilar products (e.g. Etacept, Intacept) in bioanalytical assays.

Our anti-etanercept antibodies can be used to develop pharmacokinetic (PK) bridging assays to measure either free drug or total drug. This antibody, in monovalent Fab format, is recommended as the capture antibody. For an assay to detect free drug, use with the inhibitory antibody, clone AbD25940\_hIgG4Pro ([HCA277P](#)) for detection; for an assay to detect total drug use with p75-Fc-fusion specific clone AbD26183\_hIgG1 ([HCA279P](#)) for detection.

Etanercept (Enbrel<sup>®</sup>) is a dimeric fusion protein consisting of the extracellular ligand-binding portion of the human 75 kDa (p75) TNFR2 linked to the Fc portion of human IgG1. The Fc component of etanercept contains the CH2 domain, the CH3 domain and hinge region, but not the CH1 domain of IgG1. It is used to treat autoimmune diseases by interfering with TNF, acting as a TNF inhibitor. It binds to both TNF $\alpha$  and TNF $\beta$  (Lymphotoxin A).

[View a summary of all anti-etanercept antibodies.](#)

**Affinity**

The intrinsic affinity of the monovalent form of Human anti Etanercept antibody, clone AbD25939 is  $K_D = 3.2$  nM as measured by real time, label free molecular interaction analysis on immobilized Etanercept.

**ELISA**

This product may be used in a direct or indirect ELISA or as capture antibody for etanercept.

Protocol: [PK bridging ELISA to measure free drug](#)

Protocol: [PK bridging ELISA to measure total drug](#)

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Guarantee**

12 months from date of despatch.

**Acknowledgements**

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.

Enbrel<sup>®</sup> is a registered trademark of Immunex Corporation.

His-tag is a trademark of EMD Biosciences.

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/HCA276">https://www.bio-rad-antibodies.com/SDS/HCA276</a> 10040
<b>Licensed Use</b>	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only.
<b>Technical Advice</b>	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a> .

## Related Products

### Recommended Useful Reagents

[HUMAN ANTI ETANERCEPT \(HCA277\)](#)  
[HUMAN ANTI ETANERCEPT \(HCA279\)](#)  
[HISPEC ASSAY DILUENT \(BUF049A\)](#)  
[HUMAN ANTI ETANERCEPT:HRP \(HCA277P\)](#)  
[HUMAN ANTI ETANERCEPT:HRP \(HCA279P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M420993:230706'

Printed on 01 Mar 2024