

Datasheet: HCA280

BATCH NUMBER 170447

Description:	HUMAN ANTI DENOSUMAB
Specificity:	DENOSUMAB
Other names:	PROLIA, XGEVA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD26295_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL phage display library and expressed in a human cell line - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Denosumab

Specificity	Human Anti-Denosumab Antibody, clone AbD26295_hlgG1 is a paratope specific anti-idiotypic antibody that specifically recognizes the free human monoclonal antibody denosumab. The antibody does not recognize free RANKL (receptor activator of nuclear factor kappa-B ligand) or denosumab in complex with human RANKL and can be used to measure free denosumab and biosimilar products. Clone AbD26295_hlgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a reference standard in an anti-drug antibody (ADA) assay. A pair of anti denosumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody, in full immunoglobulin format, is recommended as the detection antibody paired with clone AbD26296 (HCA288) in monovalent Fab format as the capture antibody. Denosumab (Prolia, Xgeva) is a fully human monoclonal antibody (IgG2/kappa) for the treatment of osteoporosis, treatment-induced bone loss, bone metastases, multiple myeloma, and giant cell tumor of bone. The drug specifically binds to human RANKL a protein that acts as the primary signal to promote bone removal/resorption and prevents RANKL from activating its receptor, RANK, on the surface of osteoclasts and their precursors. View a summary of all anti-denosumab antibodies.
Affinity	The intrinsic affinity of the monovalent form of this antibody is $K_D=1.8$ nM as measured by real time, label-free molecular interaction analysis on immobilized Denosumab.
ELISA	This product may be used in direct or indirect ELISA and in an ELISA bridging assay as a detection antibody together, when conjugated to HRP, with HCA288 as the capture reagent. It is also suitable as a fully human control or calibrator antibody in an ADA bridging ELISA. Protocol: PK bridging ELISA to measure free drug. Protocol: ADA bridging ELISA.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. Prolia and Xgeva are trademarks of Amgen Inc.
Health And Safety	Material Safety Datasheet documentation #10040 available at:

Information <https://www.bio-rad-antibodies.com/SDS/HCA280>

Licensed Use For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[HUMAN ANTI DENOSUMAB \(HCA281\)](#)

[HUMAN ANTI DENOSUMAB \(HCA282\)](#)

[HUMAN ANTI DENOSUMAB \(HCA288\)](#)

[HUMAN ANTI DENOSUMAB:HRP \(HCA280P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M430089:240503'

Printed on 15 Jul 2025