

Datasheet: HCA282

BATCH NUMBER 165306

Description:	HUMAN ANTI DENOSUMAB
Specificity:	DENOSUMAB
Other names:	PROLIA, XGEVA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD26862_hlgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Denosumab

Specificity	<p>Human Anti-Denosumab Antibody, clone AbD26862_hlgG1 is a paratope specific anti-idiotypic antibody that specifically recognizes the free human monoclonal antibody denosumab. The antibody does not recognize free RANKL (receptor activator of nuclear factor kappa-B ligand) or denosumab in complex with human RANKL and can be used to measure free denosumab and biosimilar products in bioanalytical assays.</p> <p>Clone AbD26862_hlgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a reference standard in an anti-drug antibody (ADA) assay.</p> <p>Denosumab (Prolia, Xgeva) is a fully human monoclonal antibody (IgG2/kappa) for the treatment of osteoporosis, treatment-induced bone loss, bone metastases, multiple myeloma, and giant cell tumor of bone. The drug specifically binds to human RANKL a protein that acts as the primary signal to promote bone removal/resorption and prevents RANKL from activating its receptor, RANK, on the surface of osteoclasts and their precursors.</p> <p>View a summary of all anti-denosumab antibodies.</p>
Affinity	The intrinsic affinity of the monovalent form of this antibody is $K_D=0.8$ as measured nM by real time, label-free molecular interaction analysis on immobilized Denosumab.
ELISA	This antibody is fully human and can be used as a reference standard in an ADA assay. Protocol: ADA bridging ELISA .
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>Prolia and Xgeva are trademarks of Amgen Inc.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA28210040</p>
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only

Technical Advice

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[HUMAN ANTI DENOSUMAB \(HCA280\)](#)

[HUMAN ANTI DENOSUMAB \(HCA281\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M417427:230328'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)