

Datasheet: HCA285

BATCH NUMBER 148386

Description:	HUMAN ANTI HUMAN IgG1 (Fc)
Specificity:	IgG1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD27686
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL PLATINUM® phage display library. Expressed in E. coli and purified using purified using NiNTA affinity chromatography. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Recombinant human IgG1 kappa

Specificity **Human Anti-Human IgG1 Fc Antibody clone AbD27686** is a recombinant monovalent Fab antibody, that recognizes human IgG1 Fc, with an affinity of 76 nM. It has been shown to have a greatly reduced affinity (89-fold) for IgG from cynomolgus macaque compared to human IgG1. It also has strongly reduced binding to rhesus macaque IgG. It does not recognize other human immunoglobulin isotypes (IgG2, IgG3, IgG4, IgA, IgE, IgM) or mouse IgG2a.

This antibody has been generated by Bio-Rad using the HuCAL® antibody library and CysDisplay® technology. It is sold for research use, with implicit permission for use for the detection of therapeutic antibodies from samples obtained from experimental animals. When tested by Bio-Rad in an antigen capture pharmacokinetic assay format, this antibody performs similarly to Mouse Anti-Human IgG (CH2 Domain) Antibody clone R10Z8E9 ([MCA5748G](#)) in 10% cynomolgus macaque serum. It recognizes a different epitope on human IgG1 Fc than clone R10Z8E9.

[View a summary of anti-human IgG Fc specific antibodies for the detection of therapeutic antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=76$ nM by real time, label free molecular interaction analysis on immobilized hlgG1 kappa
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA28510094
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[MOUSE ANTI PENTA HISTIDINE TAG:HRP \(MCA5995P\)](#)
[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371723:200612'

Printed on 08 Nov 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)