

Datasheet: HCA285

Description:	HUMAN ANTI HUMAN IgG1 (Fc)
Specificity:	IgG1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD27686
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL PLATINUM phage display library. Expressed in <i>E. coli</i> and purified using purified using NiNTA affinity chromatography. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.
Source	<i>E.coli</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Recombinant human IgG1 kappa

Specificity	<p>Human Anti-Human IgG1 Fc Antibody clone AbD27686 is a recombinant monovalent Fab antibody, that recognizes human IgG1 Fc, with an affinity of 76 nM. It has been shown to have a greatly reduced affinity (89-fold) for IgG from cynomolgus macaque compared to human IgG1. It also has strongly reduced binding to rhesus macaque IgG. It does not recognize other human immunoglobulin isotypes (IgG2, IgG3, IgG4, IgA, IgE, IgM) or mouse IgG2a.</p> <p>This antibody has been generated by Bio-Rad using the HuCAL® antibody library and CysDisplay® technology. It is sold for research use, with implicit permission for use for the detection of therapeutic antibodies from samples obtained from experimental animals. When tested by Bio-Rad in an antigen capture pharmacokinetic assay format, this antibody performs similarly to Mouse Anti-Human IgG (CH2 Domain) Antibody clone R10Z8E9 (MCA5748G) in 10% cynomolgus macaque serum. It recognizes a different epitope on human IgG1 Fc than clone R10Z8E9.</p> <p>View a summary of anti-human IgG Fc specific antibodies for the detection of therapeutic antibodies.</p>
Affinity	The intrinsic affinity of the monovalent form of this antibody is $K_D=76$ nM as measured by real time, label free molecular interaction analysis on immobilized hIgG1 kappa.
References	1. Thisted, T. <i>et al.</i> (2024) VISTA checkpoint inhibition by pH-selective antibody SNS-101 with optimized safety and pharmacokinetic profiles enhances PD-1 response. Nat Commun. 15 (1): 2917.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>His-tag is a trademark of EMD Biosciences.</p>
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA28510040
Licensed Use	For <i>in vitro</i> research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody

technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[MOUSE ANTI PENTA HISTIDINE TAG:HRP \(MCA5995P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M428451:240301'

Printed on 08 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)