

Datasheet: HCA285P

BATCH NUMBER 165952

Description:	HUMAN ANTI HUMAN IgG1 (Fc):HRP
Specificity:	IgG1
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AbD27686
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/100 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL PLATINUM phage display library. Expressed in <i>E. coli</i> and purified using purified using NiNTA affinity chromatography. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. Conjugated to horseradish peroxidase (HRP) - liquid.
Source	E.coli
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.1 mg/ml

Immunogen	Recombinant human IgG1 kappa
Specificity	<p>Human Anti-Human IgG1 Fc Antibody clone AbD27686 is a recombinant monovalent Fab antibody, that recognizes human IgG1 Fc, with an affinity of 76 nM. It has been shown to have a greatly reduced affinity (89-fold) for IgG from cynomolgus macaque compared to human IgG1. It also has strongly reduced binding to rhesus macaque IgG. It does not recognize other human immunoglobulin isotypes (IgG2, IgG3, IgG4, IgA, IgE, IgM) or mouse IgG2a.</p> <p>This antibody has been generated by Bio-Rad using the HuCAL® antibody library and CysDisplay® technology. It is sold for research use, with implicit permission for use for the detection of therapeutic antibodies from samples obtained from experimental animals. When tested by Bio-Rad in an antigen capture pharmacokinetic assay format, this antibody performs similarly to Mouse Anti-Human IgG (CH2 Domain) Antibody clone R10Z8E9 (MCA5748G) in 10% cynomolgus macaque serum. It recognizes a different epitope on human IgG1 Fc than clone R10Z8E9.</p> <p>View a summary of anti-human IgG Fc specific antibodies for the detection of therapeutic antibodies.</p>
Affinity	The intrinsic affinity of the monovalent form of this antibody is $K_D=76$ nM as measured by real time, label free molecular interaction analysis on immobilized hIgG1 kappa.
Storage	<p>This product is shipped frozen.</p> <p>When ready to use, thaw and aliquot the sample as needed. Store aliquots at -70°C, if available, otherwise store at -20°C. It is not recommended to keep aliquots at 4°C for more than one week.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA285P10094</p>
Licensed Use	For <i>in vitro</i> research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Useful Reagents

HISPEC ASSAY DILUENT (BUF049A)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M416810:230302'

Printed on 08 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)