

Datasheet: HCA288

BATCH NUMBER 163168

Description:	HUMAN ANTI DENOSUMAB
Specificity:	DENOSUMAB
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD26296
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL® phage display library, expressed in *E. coli*. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Denosumab

Specificity

Human Anti-Denosumab Antibody, clone AbD26296 is a paratope specific,

anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug denosumab. The antibody does not recognize free RANKL (receptor activator of nuclear factor kappa-B ligand), or denosumab in complex with human RANKL and can be used to measure the levels of free denosumab and biosimilar products in bioanalytical assays.

A pair of anti-denosumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody, in monovalent Fab format, is recommended as the capture antibody, paired with an anti-denosumab antibody in full immunoglobulin format, clone AbD26295_hIgG1 ([HCA280](#)) as the detection antibody.

Denosumab (Prolia, Xgeva) is a fully human monoclonal antibody (IgG2/kappa) for the treatment of osteoporosis, treatment-induced bone loss, bone metastases, multiple myeloma, and giant cell tumor of bone. The drug specifically binds to human RANKL, a protein that acts as the primary signal to promote bone removal/resorption and prevents RANKL from activating its receptor, RANK, on the surface of osteoclasts and their precursors.

[View a summary of all anti-denosumab antibodies.](#)

Affinity	The monovalent intrinsic affinity of AbD26296 was measured as $K_D = 7.1$ nM by real time, label free molecular interaction analysis on immobilized denosumab
ELISA	Clone AbD26296 can be used in a direct or indirect ELISA system or as capture antibody for denosumab in a bridging ELISA together with HCA280P (AbD26295_hIgG1) as the detection reagent. Protocol: PK bridging ELISA to measure total drug
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.</p> <p>Prolia and Xgeva are registered trademarks of Amgen Inc.</p> <p>His-tag is a registered trademark of EMD Biosciences.</p>
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA288 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any

form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory	For research purposes only
-------------------	----------------------------

Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual
-------------------------	---

Related Products

Recommended Useful Reagents

[MOUSE ANTI PENTA HISTIDINE TAG:HRP \(MCA5995P\)](#)

[HUMAN ANTI DENOSUMAB \(HCA281\)](#)

[HUMAN ANTI DENOSUMAB \(HCA282\)](#)

[HUMAN ANTI DENOSUMAB \(HCA280\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI DENOSUMAB:HRP \(HCA280P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391801:211015'

Printed on 01 Mar 2024