

Datasheet: HCA290

Description:	HUMAN ANTI GOLIMUMAB
Specificity:	GOLIMUMAB
Other names:	SIMPONI
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD25451_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Golimumab

Specificity

Human Anti-Golimumab Antibody, clone AbD25451_hIgG1 is an anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug golimumab. This antibody binds to both free golimumab and to golimumab bound to its target tumor necrosis factor alpha (TNF α). Clone AbD25451_hIgG1 does not inhibit the binding of golimumab to its target TNF α and can be used to measure total golimumab and biosimilar products in bioanalytical assays.

Clone AbD25451_hIgG1 is an affinity-matured variant of AbD20897 ([HCA242](#)), with about 17-fold affinity improvement.

A pair of anti golimumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure total drug. This antibody, in full immunoglobulin format, is recommended as the detection antibody paired with clone AbD25455 ([HCA289](#)) as the capture antibody.

Golimumab, marketed under the brand name Simponi is a recombinant human IgG1/kappa monoclonal antibody approved for the treatment of rheumatoid arthritis, psoriatic arthritis and ankylosing spondylitis. This therapeutic antibody, directed against TNF α , acts by blocking the binding of TNF α to its receptors, resulting in a down-regulation of the inflammatory response associated with autoimmune diseases.

[View a summary of all anti-golimumab antibodies](#)

Affinity	The monovalent intrinsic affinity of AbD25451_hlgG1 was measured as $K_D = 0.4$ nM by real time, label free molecular interaction analysis on immobilized golimumab
ELISA	This product may be used in a direct or indirect ELISA or as detection antibody for golimumab in a bridging ELISA together with HCA289 (AbD25455, Fab-FH) as the capture reagent. Protocol: PK bridging ELISA to measure total drug
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Simponi is a registered trademark of Janssen Biotech, Inc. in the USA and Merck & Co, Inc. in Europe
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: 10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HUMAN ANTI GOLIMUMAB \(HCA289\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com
'M335247:181204'

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

Printed on 11 Oct 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)