

## Datasheet: HCA298

**BATCH NUMBER 171095**

<b>Description:</b>	HUMAN ANTI PEMBROLIZUMAB
<b>Specificity:</b>	PEMBROLIZUMAB
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD30689_hlgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A.
<b>Source</b>	HKB-11
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml.
<b>Immunogen</b>	Pembrolizumab
<b>Specificity</b>	<b>Human Anti-Pembrolizumab Antibody, clone AbD30689_hlgG1</b> is a paratope specific,

inhibitory anti-idiotypic antibody that specifically recognizes the free humanized monoclonal antibody pembrolizumab. The antibody does not recognize recombinant human programmed cell death 1 (PD-1) or pembrolizumab in complex with recombinant human PD-1.

Clone AbD30689\_hIgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a reference standard in an anti-drug antibody (ADA) assay.

Pembrolizumab (Keytruda) is a humanized (IgG4/kappa) antibody with the heavy chain mutation S228P (IgG4-Pro). Pembrolizumab is used to treat inoperable or metastatic melanoma, metastatic non-small cell lung cancer (NSCLC), and as a second-line treatment for head and neck squamous cell carcinoma (HNSCC). For NSCLC, pembrolizumab is a first line treatment if the cancer overexpresses PDL1 and has no mutations in EGFR or ALK. Pembrolizumab is a checkpoint inhibitor and acts as an immunomodulator by blocking ligand activation of PD-1 receptor on T cells, thereby activating T cells to attack the cancer. It blocks PD-1 interaction with PD-1 ligand (PD-L1) and PD-1 ligand 2 (PD-L2).

[View a summary of all anti-pembrolizumab antibodies](#)

---

**Affinity** The intrinsic affinity of the monovalent form of Human anti Pembrolizumab antibody, clone AbD30689\_hIgG1 is  $K_D = 0.6$  nM as measured by real time, label free molecular interaction analysis on immobilized Pembrolizumab.

---

**ELISA** Human anti Pembrolizumab antibody, clone AbD30689\_hIgG1 can be used in a direct or indirect ELISA or as a detection antibody for pembrolizumab. It may also be used as a surrogate positive control to generate a standard curve for an anti-drug antibody assay.

Protocol: [anti-drug antibody assay](#).

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at  $-20^{\circ}\text{C}$  on receipt. When thawed, aliquot the sample as needed. Keep aliquots at  $2-8^{\circ}\text{C}$  for short term use (up to 4 weeks) and store the remaining aliquots at  $-20^{\circ}\text{C}$ .

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch.

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details. Keytruda is a trademark of Merck & Co.

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA298>

---

**Licensed Use** For *in vitro* research purposes and for commercial applications for the provision of *in vitro*

testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Regulatory** For research purposes only.

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

---

## Related Products

### Recommended Negative Controls

[RECOMBINANT HUMAN IgG1 LAMBDA ALLOTYPED G1m3 \(HCA049\)](#)

### Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[HUMAN ANTI PEMBROLIZUMAB \(HCA296\)](#)

[HUMAN ANTI PEMBROLIZUMAB \(HCA297\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M430095:240503'

**Printed on 21 Nov 2025**