

## Datasheet: HCA300

**BATCH NUMBER 1804**

<b>Description:</b>	HUMAN ANTI NIVOLUMAB
<b>Specificity:</b>	NIVOLUMAB
<b>Other names:</b>	Opdivo
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD30255_hIgG1
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (lambda light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
<b>Immunogen</b>	Nivolumab
<b>Specificity</b>	<b>Human Anti-Nivolumab Antibody, clone AbD30255_hIgG1</b> is a paratope specific,

inhibitory anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug nivolumab. The antibody does not recognize recombinant human programmed cell death 1 (PD-1) or nivolumab in complex with recombinant human PD-1 and can be used to measure free nivolumab levels in serum from patients.

Clone AbD30255\_hIgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a reference standard in an anti-drug antibody (ADA) assay.

Nivolumab (Opdivo) is a fully human antibody (IgG4/kappa), with the heavy chain mutation S228P (IgG4-Pro). It is used as a first line treatment for inoperable or metastatic melanoma in combination with ipilimumab if the cancer does not have a mutation in BRAF. It is also approved for treatment of advanced non-small cell lung cancer, advanced renal cell carcinoma, classical Hodgkin's lymphoma and advanced head and neck squamous cell carcinoma. Nivolumab is a checkpoint inhibitor and acts as an immuno-modulator by blocking ligand activation of PD-1 receptor on T cells, thereby activating T cells to attack the cancer. Nivolumab blocks PD-1 interaction with PD-1 ligand (PD-L1) and PD-1 ligand 2 (PD-L2).

[View a summary of all anti-nivolumab antibodies](#)

---

<b>Affinity</b>	The monovalent intrinsic affinity of AbD30255_hIgG1 was measured as $K_D = 0.5$ nM by real time, label free molecular interaction analysis on immobilized nivolumab.
<b>ELISA</b>	Clone AbD30255_hIgG1 can be used in a direct or indirect ELISA. Furthermore, it may be used as a calibrator for the set-up of an anti-drug antibody assay.  <a href="#">Protocol: anti-drug antibody assay</a>
<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Opdivo is a trademark of Bristol-Myers Squibb Company.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10094 available at: 10094: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf</a>
<b>Licensed Use</b>	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
<b>Regulatory</b>	For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

---

## Related Products

### Recommended Negative Controls

[RECOMBINANT HUMAN IgG1 LAMBDA ALLOTYPIC G1m3 \(HCA049\)](#)

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[HUMAN ANTI NIVOLUMAB \(HCA299\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M371742:200612'

Printed on 02 Mar 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)