

Datasheet: HCA310D550

BATCH NUMBER 147935

Description:	HUMAN ANTI MOUSE IgG2a:DyLight®550
Specificity:	IgG2a
Format:	DyLight®550
Product Type:	Monoclonal Antibody
Clone:	AbD24124
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/500 - 1/10,000
Western Blotting	▪			1/1000 - 1/10,000
Immunofluorescence	▪			1/100 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Product Form

A bivalent human recombinant Fab selected from the HuCAL® phage display library, expressed in E. coli. This Fab fragment is dimerized via an inactive variant of the bacterial alkaline phosphatase (BAP) fusion protein. The antibody is tagged with DYKDDDDK and His (HHHHHH) tags. This antibody is supplied conjugated to DyLight®550 - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	DyLight®550	562	576

Preparation

Prepared by affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations	Antibody concentration 1.0 mg/ml
Specificity	Human anti Mouse IgG2a monoclonal antibody , clone AbD24124, recognizes the IgG2a isotype of the primary antibody, with no species cross reactivity against mouse IgG1 and mouse IgG2b primary antibodies. This antibody allows for more than one unlabelled mouse monoclonal antibody to be used at any one time permitting multiplexing.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. His-tag is a registered trademark of EMD Biosciences. DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA310D550 10040
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M371759:200612'

Printed on 01 Mar 2024