

Datasheet: HCA313

Description:	HUMAN ANTI ECULIZUMAB
Specificity:	ECULIZUMAB
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD32107_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.

Preparation

Purified IgG prepared by affinity chromatography on Protein A

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml

Immunogen

Eculizumab

Specificity

Human Anti-Eculizumab Antibody, clone AbD32107_hIgG1 is a paratope specific, inhibitory anti-idiotypic antibody that specifically recognizes the free humanized monoclonal antibody eculizumab. It does not recognize free human complement C5 or eculizumab in complex with C5. The antibody can be used to measure free eculizumab or eculizumab biosimilars in patient samples.

Clone AbD32107_hIgG1 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. A pair of anti-eculizumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to

measure free drug. This antibody is recommended as the detection reagent, paired with Anti-Eculizumab Antibody clone AbD32311 ([HCA312](#)) in monovalent Fab format as the capture antibody.

Eculizumab (Soliris) is a humanized IgG2/4 kappa antibody. The heavy chains are comprised of human IgG2 sequences in CH1, the hinge and the adjacent portion CH2, and human IgG4 sequences in the remaining part of CH2 and CH3. The drug is used to treat paroxysmal nocturnal hemoglobinuria (PNH), a rare acquired life-threatening disease of the blood characterized by destruction of red blood cells by the complement system. It also has approval as an orphan drug for treatment of atypical hemolytic uremic syndrome (aHUS). Eculizumab specifically binds to the terminal complement component 5 (C5) and inhibits the cleavage of C5 to C5a and C5b by the C5 convertase, which prevents the generation of the terminal complement complex C5b-9. Both C5a and C5b-9 cause the terminal complement-mediated events that are characteristic of PNH and aHUS. A 1:1 molar ratio of eculizumab to C5 is sufficient for C5 inhibition.

[View a summary of all anti-eculizumab antibodies](#)

Affinity	The monovalent intrinsic affinity of AbD32107_hlgG1 was measured as $K_D = 0.6$ nM by real time, label free molecular interaction analysis on immobilized eculizumab.
ELISA	Clone AbD32107_hlgG1 can be used in an indirect ELISA or as detection antibody for eculizumab in a bridging ELISA together with HCA312 as the capture reagent. It can also be used as a surrogate positive control or calibrator in an anti-drug antibody assay. Protocol: PK bridging ELISA to measure free drug Protocol: Anti-Drug Antibody assay
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Soliris is a trademark of Alexion Pharmaceuticals, Inc.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: 10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HUMAN ANTI ECULIZUMAB \(HCA312\)](#)

[HUMAN ANTI ECULIZUMAB \(HCA314\)](#)

[HUMAN ANTI ECULIZUMAB \(HCA315\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M371769:200612'

Printed on 11 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)