

Datasheet: HCA316

BATCH NUMBER 173130

Description:	HUMAN ANTI RANIBIZUMAB
Specificity:	RANIBIZUMAB
Other names:	Lucentis
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD29858_hIgG1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Ranibizumab

Specificity **Human Anti-Ranibizumab Antibody, clone AbD29858_hIgG1** is a non-inhibitory, anti-idiotypic antibody in full immunoglobulin format that specifically recognizes the monoclonal antibody drug ranibizumab and ranibizumab when in complex with recombinant human vascular endothelial growth factor A (VEGF-A). It does not recognize free human VEGF-A or bevacizumab. This antibody can be used to measure levels of free ranibizumab and biosimilars in patient samples.

Clone AbD29858_hIgG1 in full immunoglobulin format can be used to develop a pharmacokinetic (PK) antigen capture format assay to measure free ranibizumab that has been captured via immobilized VEGF-A.

Ranibizumab (Lucentis) is a humanized Fab fragment antibody (VH3-23/kappa1), originating from the same parental mouse antibody as bevacizumab (Avastin); it has 4 amino acid exchanges in the heavy chain and 1 amino acid exchange in the light chain compared to bevacizumab. Ranibizumab was first approved to treat wet age-related macular degeneration (AMD), a common form of age-related vision loss. It is also approved for the treatment of macular edema following retinal vein occlusion, diabetic macular edema, and diabetic retinopathy. Ranibizumab inhibits angiogenesis by binding to VEGF-A, blocking its interaction with VEGF receptors on the surface of endothelial cells, and thus preventing the subsequent growth of new blood vessels.

[View a summary of all anti-ranibizumab antibodies](#)

Affinity The intrinsic affinity of the monovalent form of Human anti Ranibizumab antibody, clone AbD29858_hIgG1 is $K_D=18$ nM as measured by real time, label-free molecular interaction analysis on immobilized VEGF/ranibizumab complex.

ELISA Human anti Ranibizumab antibody, clone AbD29858_hIgG1 can be used in a direct or indirect ELISA on Ranibizumab or as a detection reagent for free Ranibizumab captured via immobilized VEGF-A.

Protocol: [PK antigen capture ELISA to measure free drug](#) .

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Lucentis is a trademark of Genentech.

Health And Safety Material Safety Datasheet documentation #10040 available at:

Information	https://www.bio-rad-antibodies.com/SDS/HCA316
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Negative Controls

[RECOMBINANT HUMAN IgG1 LAMBDA ALLOTYPE G1m3 \(HCA049\)](#)

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[RECOMBINANT HUMAN VEGF-A \(PHP293\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M430128:240503'

Printed on 19 Mar 2026