

## Datasheet: HCA322

<b>Description:</b>	HUMAN ANTI BrdU
<b>Specificity:</b>	BrdU
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD33758kd
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/200 - 1/2000
Immunohistology - Paraffin	▪			

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation program. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Chemical

### Product Form

Recombinant chimeric human/rat IgG2a antibody selected from the HuCAL phage display library and expressed in a human cell line. In this antibody, the VH and VL regions are human while the CH1-CH3 and C-lambda regions are derived from rat IgG2a. This antibody is supplied as a liquid.

### Preparation

Purified IgG prepared by affinity chromatography on Protein A

### Source

HKB-11

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.01% Thiomersal

### Approx. Protein

Antibody concentration 0.5 mg/ml

## Concentrations

---

**Immunogen** Bromodeoxyuridine

---

**Specificity** **Human anti BrdU antibody, clone AbD33758kd** recognizes bromodeoxyuridine (known as BrdU or BrdUrd). BrdU is a synthetic thymidine analogue, which is incorporated into the newly synthesized DNA during replication instead of thymidine.

Antibody detection of incorporated BrdU in cellular DNA is extensively referenced as an accurate method to monitor cell proliferation *in vivo* and *in vitro*. In cell proliferation assays BrdU staining is coupled with the use of a dye that binds total DNA such as propidium iodide (PI). BrdU can be administered diluted in the culture medium or, *in vivo* via intraperitoneal injection, subcutaneous osmotic pump implants ([Tesfaigzi et al. 2004](#)) or in drinking water ([Moser et al. 2004](#)).

BrdU can be used as a thymidine analogue in a wide range of organisms ranging from mammalian cells, through reptiles and amphibians to invertebrate species and plants.

Human anti BrdU antibody, clone AbD33758kd, has been tested in flow cytometry. Anti rat IgG2a specific antibodies are recommended as secondary detection reagents.

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.

---

**Health And Safety Information** Material Safety Datasheet documentation #10094 available at: <https://www.bio-rad-antibodies.com/SDS/HCA322>  
10094

---

**Licensed Use** For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

---

## Related Products

## Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...) [DyLight®405](#), [DyLight®488](#), [DyLight®800](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M416865:230302'

**Printed on 25 Mar 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)