

Datasheet: HCA330

Description:	HUMAN ANTI IPILIMUMAB
Specificity:	IPILIMUMAB
Other names:	Yervoy
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD34433
Isotype:	HuCAL Fab monovalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (kappa light chain) selected from the HuCAL phage display library, expressed in *E. coli*. The antibody is tagged with a DYKDDDDK tag and a double extended Strep-tag (SAWSHPQFEKGGGSGGGSGGSAWSHPQFEK). This antibody is supplied liquid.

Preparation

StrepTactin affinity chromatography.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.01% Thiomersal.

Approx. Protein Concentrations

Antibody concentration 0.5 mg/ml.

Immunogen

Ipilimumab

Specificity

Human Anti-Ipilimumab Antibody, clone AbD34433 is a paratope specific, inhibitory anti-idiotypic antibody that specifically recognizes the monoclonal antibody ipilimumab. It does not recognize free human CTLA-4 (cytotoxic T-lymphocyte antigen-4) or ipilimumab in complex with CTLA-4. The antibody can be used to measure free ipilimumab in patient samples.

A pair of anti-ipilimumab antibodies can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. This antibody, in monovalent Fab format, is recommended as the capture antibody, paired with an HRP conjugated Anti-Ipilimumab Antibody in full immunoglobulin format, clone AbD34429ia ([HCA329P](#)) as the detection antibody.

Ipilimumab (Yervoy) is a human IgG1/kappa antibody developed from a transgenic mouse. It has been approved by the FDA for the treatment of metastatic melanoma, renal cell carcinoma, and in combination with nivolumab for the treatment of previously treated microsatellite instability-high/deficient mismatch repair (MSI-H/dMMR) metastatic colorectal cancer. Ipilimumab activates the immune system by targeting CTLA-4, a protein receptor that downregulates the immune system. The action of cytotoxic T lymphocytes (CTLs) to recognize and destroy cancer cells is subject to an inhibitory mechanism that interrupts this destruction. Ipilimumab turns off this inhibitory mechanism and allows CTLs to function.

[View a summary of all anti-ipilimumab antibodies](#)

Affinity

The intrinsic affinity of the monovalent form of Human anti Ipilimumab antibody, clone AbD34433 is $K_D = 0.3$ nM as measured by real time, label free molecular interaction analysis on immobilized Ipilimumab.

ELISA

Human anti Ipilimumab antibody, clone AbD34433 can be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with [HCA329P](#) as the detection reagent.

Protocol: [PK bridging ELISA](#).

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch.

Acknowledgements

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Yervoy is a trademark of Bristol-Myers Squibb Company.

Health And Safety

Material Safety Datasheet documentation #10094 available at:

Information	10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only.
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual .

Related Products

Recommended Secondary Antibodies

Rat Anti Synthetic Peptide DYKDDDDK TAG (MCA4764...) [HRP](#)

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI IPILIMUMAB \(HCA327\)](#)

[HUMAN ANTI IPILIMUMAB \(HCA328\)](#)

[HUMAN ANTI IPILIMUMAB \(HCA329\)](#)

[HUMAN ANTI IPILIMUMAB \(DRUG/TARGET COMPLEX\) \(HCA331\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M406775:220930'

Printed on 30 Sep 2022