

Datasheet: HCA335

BATCH NUMBER 161591

Description:	HUMAN ANTI ABATACEPT
Specificity:	ABATACEPT
Other names:	Orencia
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD37058ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Abatacept
Specificity	Human Anti-Abatacept Antibody, clone AbD37058ia is a recombinant antibody that

specifically recognizes the fusion protein drug abatacept as well as cytotoxic T-lymphocyte antigen-4 (CTLA-4), also known as CD152. The antibody binds to free abatacept and can be used to measure free abatacept or abatacept biosimilars in patient samples.

Clone AbD37058ia is a fully human recombinant monoclonal antibody with an IgG1 isotype and is suitable as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. Additionally, the antibody can be used to develop a pharmacokinetic (PK) bridging assay to measure free drug. It is recommended as the detection antibody, paired with Mouse Anti-Human CD152 Antibody, clone BNI3 ([MCA1724](#)) as the capture reagent.

Abatacept (Orencia) is a recombinant fusion protein consisting of the extra-cellular domain of human CTLA-4, fused to the modified Fc (hinge, CH2, and CH3 domains) portion of human IgG1. It is used to treat several autoimmune diseases, including rheumatoid arthritis, psoriatic arthritis and juvenile idiopathic arthritis, by inhibiting optimal T-cell activation. It binds to CD80 and CD86 on antigen-presenting cells and blocks the costimulatory interaction with CD28 on T cells.

[View a summary of all anti-abatacept antibodies](#)

Affinity	The monovalent intrinsic affinity of AbD37058ia was measured as $K_D = 1$ nM by real time, label free molecular interaction analysis on immobilized abatacept.
-----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------

ELISA	Clone AbD37058ia can be used in an indirect ELISA or as detection antibody for abatacept in a bridging ELISA together with MCA1724 as the capture reagent. It can also be used as a calibrator in an anti-drug antibody assay.
--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Protocol: [PK bridging ELISA](#)

Protocol: [ADA bridging ELISA](#)

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. Orencia is a trademark of Bristol-Myers Squibb Company.
-------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA335 10094
--------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Licensed Use For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI ABATACEPT \(HCA336\)](#)

ELISA Matched Pair - Capture Antibody

[MOUSE ANTI HUMAN CD152 \(MCA1724\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391864:211015'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)