

Datasheet: HCA344

BATCH NUMBER 172968

Description:	HUMAN ANTI PERTUZUMAB
Specificity:	PERTUZUMAB
Other names:	Perjeta
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD35668ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A.
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Pertuzumab

Specificity	<p>Human Anti-Pertuzumab Antibody, clone AbD35668ia is a paratope specific, inhibitory anti-idiotypic antibody that specifically recognizes the humanized monoclonal antibody pertuzumab. It does not recognize free human epidermal growth factor receptor 2 (HER2).</p> <p>Clone AbD35668ia is a fully human recombinant monoclonal antibody with IgG1 isotype and can be used as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 14 nM, and is one of two antibodies with different affinities offered for ADA assays.</p> <p>Pertuzumab (Perjeta) is an IgG1/kappa antibody humanized from mouse. It has been approved for the treatment of patients with metastatic HER2-positive breast cancer in combination with trastuzumab. Pertuzumab binds to the extracellular domain II of HER2 and inhibits HER2–HER3 dimerization, reducing signaling through intracellular pathways (Capelan et al. 2013).</p> <p>View a summary of all anti-pertuzumab antibodies</p>
Affinity	<p>The intrinsic affinity of the monovalent form of Human anti Pertuzumab antibody, clone AbD35668ia is $K_D = 14$ nM as measured by real time, label free molecular interaction analysis on immobilized Pertuzumab.</p>
ELISA	<p>Human anti Pertuzumab antibody, clone AbD35668ia can be used in a direct or an indirect ELISA. It can also be used as a calibrator in an anti-drug antibody assay.</p> <p>Protocol: ADA bridging ELISA.</p>
Further Reading	<p>1. Capelan, M. <i>et al.</i> (2013) Pertuzumab: new hope for patients with HER2-positive breast cancer. Ann Oncol. 24 (2): 273-282.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	<p>12 months from date of despatch</p>
Acknowledgements	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. Perjeta is a trademark of Genentech, Inc.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA344</p>
Licensed Use	<p>For <i>in vitro</i> research purposes and for commercial applications for the provision of <i>in vitro</i> testing services to support preclinical and clinical drug development. Any re-sale in any</p>

form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA341\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA342\)](#)

[HUMAN ANTI PERTUZUMAB \(HCA345\)](#)

[HUMAN ANTI PERTUZUMAB \(TZA025\)](#)

[HUMAN ANTI PERTUZUMAB:HRP \(TZA025P\)](#)

[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\) \(TZA035\)](#)

[HUMAN ANTI PERTUZUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA035P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M442698:250528'

Printed on 29 Jan 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)