

Datasheet: HCA356

BATCH NUMBER 172957

Description:	ANTI PARP1-S499-ADP-RIBOSE
Specificity:	PARP1-S499-ADP-RIBOSE
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD34251
Isotype:	IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Human/Rabbit IgG chimera (rabbit CH2 and CH3) antibody selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified recombinant IgG prepared by affinity chromatography on Protein A from a mammalian cell line
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)

Carrier Free	Yes
Approx. Protein Concentrations	Total protein concentration 0.5 mg/ml
Immunogen	APRGKS(ADPr)GAALSKKSKGQVGGK
Specificity	<p>Anti-PARP1-S499-ADP-ribose antibody, clone AbD34251, recognizes PARP1, ADP-ribosylated at Serine 499. Poly(ADP-ribose) polymerase 1 (PARP1) is an early responder to DNA damage in human cells. Upon binding to genomic lesions PARP1 is able to transfer a mono- or poly-ADP-ribose residue from nicotinamide adenine dinucleotide (NAD⁺) to their substrate proteins, for chromatin decompaction and repair factor recruitment (Bütepage et al. 2015, Suskiewicz et al. 2020).</p> <p>During the DNA damage response, serine is the primary target for PARP1 ADP-ribosylation which uses Histone PARylation factor 1 (HPF1) as an accessory factor to switch the amino-acid specificity of PARP1 from aspartate/glutamate to serine residues (Bonfiglio et al. 2017, Suskiewicz et al. 2020). PARP1/HPF1 during DNA damage response, primarily target PARP1 itself and histone H3, ribosylating these substrate proteins.</p> <p>This PARP1-S499ADPr antibody, clone AbD34251 was generated using an ADP-ribosylated PARP1S499 antigen. Characterization has demonstrated the antibody does not bind any other Ser-ADP-ribosylated peptides tested and binds via three residues before and one after the modified Serine. This clone does not detect poly-ADP-ribose in a site-specific manner (Bonfiglio et al. 2020).</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/HCA356
Regulatory	For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

Printed on 09 Mar 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)