

## Datasheet: HCA357

**BATCH NUMBER 173129**

<b>Description:</b>	ANTI H3-S10-ADP-RIBOSE
<b>Specificity:</b>	H3-S10-ADP-RIBOSE
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD33644
<b>Isotype:</b>	HuCAL Fab bivalent
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Human

### Product Form

A bivalent human recombinant Fab (lambda light chain) selected from the HuCAL phage display library. Expressed in *E. coli*. This Fab fragment is bivalent by dimerization of the bacterial alkaline phosphatase fusion protein. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as liquid.

### Preparation

Metal chelate affinity chromatography

### Source

*E. coli*

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	Total protein concentration 0.5 mg/ml
<b>Immunogen</b>	ARTKQTARKS(ADPr)TGGKAC
<b>Specificity</b>	<p><b>Anti-H3-S10-ADP-ribose antibody, clone AbD33644</b>, recognizes Histone H3, ADP-ribosylated at Serine 10. This antibody also recognizes the less abundant H3 site H3-S28-ADP-ribose, however this is with a lower affinity and is due to the similarities between sequences flanking the two H3 target sites. This clone does not detect poly-ADP-ribose in a site-specific manner (<a href="#">Bonfiglio et al. 2020</a>).</p> <p>Poly(ADP-ribose) polymerase 1 (PARP1) is an early responder to DNA damage in human cells. Upon binding to genomic lesions PARP1 is able to transfer a mono- or poly-ADP-ribose residue from nicotinamide adenine dinucleotide (NAD<sup>+</sup>) to their substrate proteins, for chromatin decompaction and repair factor recruitment (<a href="#">Bütepage et al. 2015</a>, <a href="#">Suskiewicz et al. 2020</a>).</p> <p>During the DNA damage response, serine is the primary target for PARP1 ADP-ribosylation which uses Histone PARylation factor 1 (HPF1) as an accessory factor to switch the amino-acid specificity of PARP1 from aspartate/glutamate to serine residues (<a href="#">Bonfiglio et al. 2017</a>, <a href="#">Suskiewicz et al. 2020</a>). During the DNA damage response PARP1 and HPF1 primarily target PARP1 itself and histone H3, ribosylating these targets.</p> <p>This clone has been used in Western blotting, in the context of the DNA damage response with Histone H3 Ser-ADPr sites recognized by this antibody in lysates from Wild Type and ARH3 KO cells treated with H<sub>2</sub>O<sub>2</sub> (<a href="#">Bonfiglio et al. 2020</a>).</p>
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	<p>This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <a href="http://bio-rad.com/en-us/trademarks">bio-rad.com/en-us/trademarks</a> for details.</p> <p>His-tag is a trademark of EMD Biosciences.</p>
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/HCA357">https://www.bio-rad-antibodies.com/SDS/HCA357</a></p>
<b>Regulatory</b>	For research purposes only

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M449518:260127'

**Printed on 09 Mar 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)