

Datasheet: HCA359

BATCH NUMBER 163201

Description:	HUMAN ANTI DURVALUMAB
Specificity:	DURVALUMAB
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD41310ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Drug
Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Durvalumab
Specificity	Human Anti-Durvalumab Antibody, clone AbD41310ia is a paratope specific, inhibitory

anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug durvalumab. The antibody does not recognize free PD-L1 nor durvalumab in complex with PD-L1.

Clone AbD41310ia is a fully human recombinant monoclonal antibody with IgG1 isotype and can be used as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 4.8 nM, and is one of two antibodies with different affinities offered for ADA assays.

[ADA bridging ELISA](#)

Durvalumab (Imfinzi) is a human IgG1/kappa antibody with allotype G1m3. The CH2 domain was engineered with three amino acid substitutions (L238F, L239E, P335S) to reduce the Fc-mediated effector functions to prevent depletion of PD-L1+ T-cells. It was first approved for previously treated patients with advanced bladder cancer. It was later approved for the treatment of unresectable stage 3 NSCLC (Non-Small Cell Lung Cancer) and ES-SCLC. Durvalumab binds to PD-L1 to block the interaction with the inhibitory T cell receptor PD-1. This blockade removes the suppression of T-cells that are specific for MHC-1-presented neoepitopes on tumor cells. Removal of this inhibitory signal leads to the release of granzymes and perforins by the activated T-cells, and finally to the destruction of the tumor cells.

[View a summary of all anti-durvalumab antibodies](#)

Affinity	The intrinsic affinity of the monovalent form of Human anti Durvalumab antibody, clone AbD41310ia is $K_D = 4.8$ nM as measured by real time, label free molecular interaction analysis on immobilized Durvalumab.
ELISA	Human anti Durvalumab antibody, clone AbD41310ia can be used as a calibrator in an anti-drug antibody assay. Protocol: ADA bridging ELISA .
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. Imfinzi is a trademark of Medimmune/AstraZeneca using Abgenix's XenoMouse technology.

Health And Safety Information Material Safety Datasheet documentation #10094 available at:
10094: <https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

Licensed Use For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI DURVALUMAB \(HCA360\)](#)

[HUMAN ANTI DURVALUMAB \(TZA004\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M412558:221114'

Printed on 02 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)