

Datasheet: HCA364

BATCH NUMBER 162455

Description:	HUMAN ANTI DARATUMUMAB
Specificity:	DARATUMUMAB
Other names:	Darzalex
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD45015ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL® phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Daratumumab
Specificity	Human Anti-Daratumumab Antibody, clone AbD45015ia is a paratope specific,

inhibitory anti-idiotypic antibody that specifically recognizes the monoclonal antibody daratumumab. It does not recognize free CD38 nor daratumumab in complex with CD38. The antibody can be used to measure free daratumumab in patient samples.

This anti-daratumumab antibody can be used as both capture and detection reagent to develop a pharmacokinetic (PK) bridging assay to measure free drug.

Clone AbD45015ia is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 0.9 nM, and is one of three antibodies with different affinities offered for ADA assays.

Daratumumab (Darzalex) is a human IgG1/kappa antibody generated using transgenic mouse technology. It has been approved for the treatment of patients with multiple myeloma. It has orphan drug designation for multiple myeloma, diffuse large B cell lymphoma, follicular lymphoma, and mantle cell lymphoma. The antibody drug targets CD38, overexpressed by multiple myeloma cells, and indirectly induces apoptosis. In addition, daratumumab modulates CD38 enzymatic activity, inhibiting the cyclase enzyme activity and stimulating the hydrolase activity ([Van de Donk *et al.* 2016](#)).

[View a summary of all anti-daratumumab antibodies](#)

Affinity	The monovalent intrinsic affinity of AbD45015ia was measured as $K_D = 0.9$ nM by real time, label free molecular interaction analysis on immobilized daratumumab.
ELISA	<p>Clone AbD45015ia can be used as capture and detection antibody for daratumumab in a bridging ELISA to measure free drug. It can also be used as a calibrator in an anti-drug antibody assay.</p> <p>Protocol: ADA bridging ELISA</p> <p>Protocol: PK bridging ELISA</p>
Further Reading	1. van, N.D.D. <i>et al.</i> (2016) Monoclonal antibodies targeting CD38 in hematological malignancies and beyond. Immunol Rev. 270 (1): 95-112.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.

Darzalex is a trademark of Janssen Biotech, Inc.

Health And Safety Information	Material Safety Datasheet documentation #10094 available at: https://www.bio-rad-antibodies.com/SDS/HCA364 10094
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Useful Reagents

[HUMAN ANTI DARATUMUMAB \(HCA365\)](#)

[HUMAN ANTI DARATUMUMAB \(HCA366\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M393943:220120'

Printed on 01 Mar 2024