

## Datasheet: HCA378

<b>Description:</b>	HUMAN ANTI EVOLOCUMAB
<b>Specificity:</b>	EVOLOCUMAB
<b>Other names:</b>	Repatha
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD40405ia
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Source</b>	HKB-11
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
<b>Immunogen</b>	Evolocumab

**Specificity** **Human Anti-Evolocumab Antibody, clone AbD40405ia**, is a paratope specific, inhibitory antibody that specifically recognizes the monoclonal antibody drug evolocumab. It does not recognize the drug target, human proprotein convertase subtilisin/kexin type 9 (PCSK9), nor evolocumab in complex with PCSK9.

Clone AbD40405ia is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 5 nM and is one of three offered for ADA assays. Antibody [HCA376](#) has affinity 5 nM and [HCA377](#) has affinity 1 nM.

Evolocumab (Repatha) is a fully human IgG2 lambda antibody developed using transgenic mouse technology. It has been approved for treating selected patients with high cholesterol, and for the prevention of heart attack and stroke. The drug prevents the PCSK9 protein binding to the low density lipoprotein receptor (LDLR), thereby increasing the number of LDLRs available and lowering LDL levels.

[View a summary of all evolocumab antibodies](#)

---

**Affinity** The intrinsic affinity of the monovalent form of Human anti Evolocumab antibody, clone AbD40405ia is  $K_D = 5$  nM as measured by real time, label free molecular interaction analysis on immobilized Evolocumab.

---

**ELISA** Human anti Evolocumab antibody, clone AbD40405ia can be used as a calibrator in an anti-drug antibody assay.  
Protocol: [ADA bridging ELISA](#).

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.  
Repatha is a trademark of Amgen.

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA378>  
10040

---

**Licensed Use** For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

---

## Related Products

### Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI EVOLOCUMAB \(HCA376\)](#)

[HUMAN ANTI EVOLOCUMAB \(HCA377\)](#)

[HUMAN ANTI EVOLOCUMAB \(TZA007\)](#)

[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

[HUMAN ANTI EVOLOCUMAB \(TZA008\)](#)

[HUMAN ANTI EVOLOCUMAB:HRP \(TZA008P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M421101:230706'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)