

Datasheet: HCA382

Description:	HUMAN ANTI ATEZOLIZUMAB
Specificity:	ATEZOLIZUMAB
Other names:	Tecentriq
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD40638ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (kappa light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Atezolizumab

Specificity

Human Anti-Atezolizumab Antibody, clone AbD40638ia is a paratope specific, inhibitory anti-idiotypic antibody (Type 1) that specifically recognizes the free monoclonal antibody drug atezolizumab. It does not recognize the drug target, Programmed Cell Death Ligand 1 (PD-L1), nor atezolizumab in complex with PD-L1. The antibody can be used to measure levels of free atezolizumab in patient samples.

Clone AbD40638ia is a fully human recombinant monoclonal antibody with IgG1 isotype and can be used as a surrogate positive control or reference standard in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 0.4 nM, and is one of three Type 1 antibodies with different affinities offered for ADA assays.

Clone AbD40638ia is also available in a monovalent Fab format with a SpyTag2 at the C-terminus of the Fab antibody heavy chain, enabling the user to couple it to a SpyCatcher reagent for conversion to alternative formats in less than one hour under product code [TZA011](#).

Atezolizumab (Tecentriq) is a human monoclonal IgG1 antibody that binds directly to Programmed Death Ligand 1 (PD-L1, CD274) and provides a dual blockade of both the PD-1 and B7.1 receptors on T-cells and antigen presenting cells. It is one of the named "checkpoint inhibitor" drugs. PD-L1/PD-1 mediated inhibition of the immune response is thereby released, including reactivation of the anti-tumour immune response without induction of antibody-dependent cellular cytotoxicity. The drug is used to treat urothelial carcinoma, non-small cell lung cancer, triple-negative breast cancer, small cell lung cancer and hepatocellular carcinoma. It is used either on its own or in combination for cancers which are advanced or metastatic.

[View a summary of all anti-atezolizumab antibodies](#)

Affinity

The intrinsic affinity of the monovalent form of Human anti Atezolizumab antibody, clone AbD40638ia is $K_D = 0.4$ nM as measured by real time, label free molecular interaction analysis on immobilized Atezolizumab.

ELISA

Human anti Atezolizumab antibody, clone AbD40638ia can be used in a direct and indirect ELISA. It can also be used as a calibrator in an anti-drug (ADA) antibody assay.

Protocol: [ADA bridging ELISA](#).

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](https://www.bio-rad.com/en-us/trademarks) for details.
Tecentriq is a trademark of Roche/ Genentech..

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA382>
10040

Licensed Use For *in vitro* research purposes and for commercial applications for the provision of *in vitro* testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI ATEZOLIZUMAB \(HCA380\)](#)

[HUMAN ANTI ATEZOLIZUMAB \(HCA381\)](#)

[HUMAN ANTI ATEZOLIZUMAB \(TZA011\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M428552:240301'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)