

Datasheet: HCA432

**BATCH NUMBER 171368**

<b>Description:</b>	HUMAN ANTI IXEKIZUMAB
<b>Specificity:</b>	IXEKIZUMAB
<b>Other names:</b>	Taltz
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD60570ia
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This product is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Source</b>	HKB-11
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
<b>Immunogen</b>	Ixekizumab

**Specificity** **Human Anti-Ixekizumab Antibody, clone AbD60570ia** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug ixekizumab, and its biosimilar products. It does not recognize the drug target, human Interleukin 17 (IL-17A), nor ixekizumab in complex with IL-17A. This antibody can be used to measure free ixekizumab levels in serum from patients.

HCA432 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 0.2nM and is one of the four antibodies offered for ADA assays. Antibody [HCA433](#) has an affinity of 0.9nM, [HCA434](#) has an affinity of 0.6nM and [HCA435](#) has an affinity of 0.6nM.

Ixekizumab (trade name Taltz) is a humanized monoclonal antibody of the IgG4/kappa isotype that binds to interleukin 17 (IL-17A), a pro-inflammatory cytokine, and blocks its action. Blocking its effect is beneficial for use in inflammatory conditions. Ixekizumab is indicated for the treatment of moderate-to-severe plaque psoriasis and as a second-line therapy for active psoriatic arthritis.

---

**Affinity** The monovalent intrinsic affinity of Human Anti-Ixekizumab Antibody, clone AbD60570ia was measured as  $K_D = 0.2 \text{ nM}$  by real time, label-free molecular interaction analysis on immobilized ixekizumab.

---

**ELISA** Clone AbD60570ia can be used in indirect ELISA to detect ixekizumab. It can also be used as a surrogate positive control or a calibrator in an anti-drug antibody (ADA) assay, together with [HCA433](#), [HCA434](#) and [HCA435](#).

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at  $-20^\circ\text{C}$  on receipt. When thawed, aliquot the sample as needed. Keep aliquots at  $2-8^\circ\text{C}$  for short term use (up to 4 weeks) and store the remaining aliquots at  $-20^\circ\text{C}$ .

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.  
Taltz is a registered trademark of Eli Lilly.

---

**Health And Safety Information** Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA432>

---

**Licensed Use** For *in vitro*. research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

---

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI IXEKIZUMAB \(HCA433\)](#)

[HUMAN ANTI IXEKIZUMAB \(HCA434\)](#)

[HUMAN ANTI IXEKIZUMAB \(HCA435\)](#)

[HUMAN ANTI IXEKIZUMAB \(TZA0105\)](#)

[HUMAN ANTI IXEKIZUMAB \(TZA0106\)](#)

[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0106P\)](#)

[HUMAN ANTI IXEKIZUMAB \(TZA0107\)](#)

[HUMAN ANTI IXEKIZUMAB:HRP \(TZA0107P\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M444935:250814'

**Printed on 21 Nov 2025**