

Datasheet: HCA443

**BATCH NUMBER 173512**

<b>Description:</b>	HUMAN ANTI ISATUXIMAB
<b>Specificity:</b>	ISATUXIMAB
<b>Other names:</b>	Sarclisa
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD64597ia
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Source</b>	HKB-11
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Antibody concentration 0.5 mg/ml
<b>Immunogen</b>	Isatuximab

**Specificity**

**Human Anti-Isatuximab Antibody, clone AbD64597ia** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug isatuximab. It does not recognize the drug target, CD38, nor isatuximab in complex with CD38. This antibody can be used to measure free isatuximab levels in serum from patients.

HCA443 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 0.4nM and is one of the three antibodies offered for ADA assays. Antibody [HCA441](#) has an affinity of 0.2nM and [HCA442](#) has an affinity of 0.2nM. Isatuximab (trade name Sarclisa) is a chimeric human/mouse IgG1/kappa, produced using recombinant DNA technology in Chinese hamster ovary (CHO) cells. Isatuximab is an IgG1-derived monoclonal antibody targeted against CD38. Its binding towards CD38 positive myeloma cells results in several downstream effects, including direct apoptosis of the affected cell and activation of immune mechanisms including antibody-dependent cell-mediated cytotoxicity (ADCC), antibody-dependent cellular phagocytosis (ADCP), and complement dependent cyto-toxicity (CDC), all of which result in potent anti-tumour activity. Isatuximab is indicated for the treatment of multiple myeloma in adults.

[View a summary of all anti-isatuximab antibodies](#)

---

**Affinity**

The monovalent intrinsic affinity of Human anti-Isatuximab Antibody, clone AbD64597ia was measured as  $K_D = 0.4 \text{ nM}$  by real time, label-free molecular interaction analysis on immobilized isatuximab.

---

**ELISA**

Clone AbD64597ia can be used in indirect ELISA to detect isatuximab. It can also be used as a surrogate positive control or a calibrator in an anti-drug antibody (ADA) assay, together with [HCA441](#) and [HCA442](#).

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at  $-20^\circ\text{C}$  on receipt. When thawed, aliquot the sample as needed. Keep aliquots at  $2-8^\circ\text{C}$  for short term use (up to 4 weeks) and store the remaining aliquots at  $-20^\circ\text{C}$ .

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See [bio-rad.com/en-us/trademarks](http://bio-rad.com/en-us/trademarks) for details.

Sarclisa is a registered trademark of ImmunoGen Inc, Sanofi, and Sanofi-Aventis U.S LLC.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA443>

**Licensed Use** For *in vitro*. research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

---

**Regulatory** For research purposes only

---

**Technical Advice** Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

---

## Related Products

### Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

[HUMAN ANTI ISATUXIMAB \(HCA441\)](#)

[HUMAN ANTI ISATUXIMAB \(HCA442\)](#)

[HUMAN ANTI ISATUXIMAB \(TZA0124\)](#)

[HUMAN ANTI ISATUXIMAB \(TZA0125\)](#)

[HUMAN ANTI ISATUXIMAB:HRP \(TZA0125P\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M450739:260519'

**Printed on 22 Jun 2026**